PCT

(33) Priority Country:

WORLD INTELLECTUAL PROPERTY ORGANIZATION International Bureau



INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification 4:

G05D 7/06

A1

(11) International Publication Number: WO 86/ 04698

(43) International Publication Date: 14 August 1986 (14.08.86)

US

(21) International Application Number: PCT/US86/00258

(22) International Filing Date: 4 February 1986 (04.02.86)

(31) Priority Application Number: 698,434

(32) Priority Date: 5 February 1985 (05.02.85)

(05) 1 1.01 kg 1705 (05.02.05

(71) Applicant: THE SCOPAS TECHNOLOGY COM-PANY, INC. [US/US]; 60 East 42nd Street, New York, NY 10165 (US).

(72) Inventors: JEFFERIS, Raymond, P., III; 276 Hillcrest Rd., Wayne, PA 19087 (US). ENGLER, Philip, V.; 21 Maple Street, Tarrytown, NY 10591 (US). ROSENBLATT, Aaron, A.; 32 West 76th Street, New York, NY 10591 (US).

(74) Agents: SINDER, Stuart, J. et al.; Kenyon and Kenyon, One Broadway, New York, NY 10004 (US).


(81) Designated States: AU, BR, DK, FI, HU, IT (European patent), JP, KR, NO.

Published

With international search report.

(54) Title: GAS STERILANT SYSTEM

(57) Abstract



A system for treating articles, preferably with a sterilizing gas. The system includes a chamber (10) into which the articles are received and valves (V2, V1) for supplying the sterilizing gas to the chamber and for removing the gas from the chamber after a predetermined time period. The sterilizing gas is generated on site from at least two components, thus minimizing problems in the transportation of the gas to the location. The sterilizing gas generated on site is preferably chlorine dioxide and the two components may be chlorine gas and sodium chlorite. The system includes a progammed microprocessor controller (100) for controlling the valves executing a predetermined sequence of instructions. The predetermined sequence of instructions define a state diagram for the system having a plurality of successive states. In order to provide for system safety, the controller preferably employs a plurality of abort states to which the system returns in the event of a failure. Depending on the nature of the failure, the system automatically moves to the proper abort state.

FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

			•		
AT	Austria	GA	Gabon ·	MR	Mauritania
ΑU	Australia	GB	United Kingdom	MW	Malawi
BB	Barbados	HU	Hungary	NL	Netherlands
BE	Belgium	П	Italy	NO	Norway
BG	Bulgaria	JP	Japan	RO.	Romania
BR	Brazil	KP	Democratic People's Republic		
CF	Central African Republic		of Korea	SD	Sudan
œ	Congo	KR		SE	Sweden
CH	Switzerland		Republic of Korea	SN	Senegal
CM	•	LI	Liechtenstein	รบ	Soviet Union
	Cameroon	LK	Sri Lanka	TD	Chad
DE	Germany, Federal Republic of	LU	Luxembourg	TG	Togo ·
DK	Denmark	MC	Monaco	us	United States of America
FI	Finland	MG	Madagascar		
FR	France	ML	Mali		

GAS STERILANT SYSTEM

BACKGROUND OF THE INVENTION

This application is related to copending applications

5 Serial Nos. 435,331 filed October 19, 1982 and 601,443, filed April 18, 1984, the disclosures of which are hereby incorporated herein by reference.

The present invention relates to systems for delivering a gas to a confined chamber and to systems for

- sterilizing substances and articles and particularly to systems using a sterilizing gas to sterilize articles, for example medical apparatus such as utensils and instruments which may have been contaminated by foreign substances. The system of the present invention can
- also be used to sterilize non-medical articles and substances, as required. The system of the present invention relates particularly to a gas sterilizing system wherein two components which react to provide sterilizing amounts of a gas are combined in the field
- 20 by the apparatus of the present invention. This allows the components which react to form the sterilizing gas to be shipped separately, which minimizes the possibility of accidents.

In particular, the present invention relates to a system 25 using chlorine dioxide as the sterilizing gas. Chlorine dioxide gas is both unstable and toxic to humans. For

example, chlorine dioxide gas, will, over time, decompose into its constituent parts and accordingly, it cannot be transported easily. It is therefore undesirable to transport chlorine dioxide gas.

5 Moreover, chlorine dioxide gas is somewhat explosive and also has a propensity to undergo catalytic decomposition. The components which react to form chlorine dioxide gas (e.g., sodium chlorite and chlorine gas), however, may be transported relatively easily and reacted on site to provide the sterilizing gas chlorine dioxide.

Prior systems have typically used ethylene oxide gas as a sterilant. For example, the castle 4040 ethylene oxide sterilizer manufactured by Sybron Corporation,

15 Medical Products Division, is an example of such a prior system. Although ethylene oxide has been used as a sterilizing gas in the prior systems, chlorine dioxide is a preferred sterilant.

Furthermore, the systems used in the past have typically been of rather simple design and have not included advanced means for maintaining the reliability of the devices and safeguarding against accidents.

Additionally, these systems have not provided a great deal of redundancy so that if a component of the system failed, manual intervention or service personnel was required to correct the failure before the sterilizing process could continue.

SUMMARY OF THE INVENTION

It is an object of the present invention to provide a sterilizing system which uses a gas having bacteriocidal, sporicidal, fungicidal and/or viricidal properties to sterilize articles.

It is a further object of the invention to provide a

sterilizing system in which at least two components which react to provide a sterilizing gas are reacted on site within the apparatus of the present invention to provide effective amounts of the sterilizing gas.

5 It is yet a further object of the present invention to provide a sterilizing system wherein the sterilant is chlorine dioxide gas.

It is still a further object of the present invention to provide a gas sterilizing system having built-in redundancy and means for maintaining the reliability and safety of the system.

It is still yet another object of the present invention to provide a gas sterilizing system which is versatile and which is controlled by a programmed microprocessor.

15 According to one embodiment of the invention, these and other objects of the present invention are achieved by a system for treating articles with a gas comprising first means for receiving a first component, second means for receiving a second component, the first and second 20 components, when reacted together, forming the gas, means for reacting the two components together for forming the gas, valve means for supplying the gas to the chamber means to treat the article in the chamber means, means for removing the gas from the chamber 25 means, electronic controller means for controlling the means for reacting, means for supplying and means for removing, comprising computer means executing a predetermined sequence of steps so as to cycle the apparatus through a series of successive states defining 30 a cycle in which the article is treated by the gas and wherein the gas is thereafter removed from the chamber means so as to render the atmosphere in the chamber

means within acceptable standards of safety.

According to another embodiment of the invention, a system for treating articles with a gas is provided comprising chamber means for receiving articles to be treated, means for supplying the gas to the chamber 5 means comprising valve means coupled to the chamber means for supplying the gas to the chamber means, means for removing the gas from the chamber means after a predetermined time interval, electronic control means receiving a plurality of electrical signals associated 10 with ones of measured parameters from the chamber means for controlling the valve means and the means for removing, the electronic control means comprising computer means for cycling the apparatus through a plurality of states in accordance with a predetermined. sequence of instructions, the computer means including means for aborting the operation of the apparatus to one of a plurality of defined failure states in response to a failure of the apparatus, the selected failure state dependent on the state in the cycle in which the failure 20 occurred.

According to still another embodiment of the invention, a system for treating articles with a gas is provided comprising chamber means for receiving articles to be treated, means for supplying the gas to the chamber 25 means comprising valve means coupled to the chamber means for supplying the gas to the chamber means, means for removing the gas from the chamber means after a predetermined time interval, electronic control means receiving a plurality of elecrical signals associated 30 with ones of measured parameters from the chamber means for controlling the valve means and the means for removing, the electronic control means comprising computer means for cycling the apparatus through a pluality of states in accordance with a predetermined 35 sequence of instructions, the computer means including memory means, and further comprising means for receiving 15

input signals from the valve means indicative of the closed or open condition of the valve means and means for transmitting output signals to the valve means to selectively open or close the valve means, image signals of the input and output signals being stored in the memory means, mask means being stored in the memory means, the computer means comparing the image signals of the input and output signals and generating an alarm signal if the input and output image signals do not agree in response to the setting of a bit in the mask means.

Other objects, features and advantages of the present invention will be apparent from the description which follows.

BRIEF DESCRIPTION OF THE DRAWINGS

The present invention will be described in greater detail with reference to the accompanying drawings in which:

- FIG. 1 is a block diagram of the overall gas sterilant 20 system according to the invention;
 - FIG. 2 is a block diagram of the sterilizing chamber and the valve and pump block of the gas sterilant system according to the present invention;
- FIG. 3 is a block diagram of the electronic control circuitry of the gas sterilant system;
 - FIG. 3A is a table of addresses used in the electronic controller of FIG. 3 and the corresponding components or signals controlled by the addresses;
- FIG. 4 is a block diagram showing how various system 30 clock frequencies and the system interrupt are derived;

- FIG. 5 is a front view of one embodiment of a control panel for the gas sterilant system showing the controller display lights and control switches;
- FIG. 6 is a state diagram for the gas sterilant system 5 according to the present invention;
 - FIG. 7 is a state output matrix corresponding to the state diagram of FIG. 6 for the gas sterilant system according to the present invention;
- FIGS. 7A and 7B are flowcharts for the sequencing 10 program for implementing the state diagram of FIG. 6;
 - FIG. 8 is a block diagram of the safety interlock arrangement for the gas sterilant system according to the present invention;
- FIG. 9 is a functional flow diagram for the software
 15 resident in the memory of the electronic controller of
 the gas sterilant system according to the present
 invention;
 - FIG. 10 is a flow diagram for timed functions of the software for the gas sterilant system;
- 20 FIG. 11 is a flow diagram for one of the timed functions of the software for the gas sterilant system;
 - FIG. 12 is a memory map of the data memory of the electronic control circuitry for the gas sterilant system according to the present invention;
- 25 FIG. 13 is a flowchart for another of the timed functions of the software of the electronic control circuitry for the gas sterilant system according to the present invention;

- FIG. 14 is a flowchart of another of the timed functions of the software for the gas sterilant system according to the present invention;
- FIG. 15 is a flowchart for another of the timed

 5 functions of the software for the gas sterilant system according to the present invention;
 - FIG. 16 is a flowchart for a program implemented in the control unit for resetting the control unit timed functions;
- 10 FIG. 17 is a flowchart for a program implemented in the control unit for reading in input data from the system according to the invention;
 - FIG. 18 is a flowchart for a program implemented in the control unit for providing a timeout alarm in the event of a component failure;
 - FIG. 19 is a flowchart for a program implemented in the control unit for providing an additional alarm in the event of a component failure;
- FIG. 20 is a flowchart for a program implemented in the control unit for writing out data to the controlled components of the system;
 - FIG. 21 is a flowchart for a program implemented in the control unit for reading in analog input data from the controlled system;
- 25 FIG. 22 is a general flowchart for a program implemented in the control unit for providing the various timed functions of the system;
 - FIG. 23 is a flowchart for part of the program of FIG.

22; and

FIG. 24 is a flowchart for a program implemented in the control unit for controlling the system outputs.

DETAILED DESCRIPTION

5 Overall System

With reference now to the drawing figures, FIG. 1 shows the overall gas sterilant system. The system comprises a sterilizing chamber 10, electronic control circuitry 100 which is preferably microprocessor controlled, valve 10 and pump block 20 and displays 110. Sensor inputs 5 including signals generated by appropriate sensors in chamber 10 and related to temperature, pressure, humidity and sterilizing gas concentration in the chamber 10 are fed from the sterilizing chamber 10 to 15 control circuitry 100. The sensor inputs include both analog signals relating to the above measured chamber parameters and certain digital signals, e.g., a signal indicative of when the temperature in the chamber has reached a desired value, to be explained in more detail 20 below. A START CYCLE switch S1 initiates operation of the system and an ABORT-RESET switch S2, as described in more detail later, is used to recycle the system states to a defined condition if an abort mode is attained, i.e., if a failure or alarm condition occurs. 25 operation of valve and pump block 20 will be described in more detail below, and includes a source of chlorine dioxide gas 22 which is produced on location from separated components, water vapor 23 and nitrogen 24. The valve and pump block is also vented to the 30 atmosphere, as shown. Valve and pump block 20 includes a number of sequenced and controlled valves and a vacuum pump for providing the necessary conditions in the sterilizing chamber at the appropriate times. of the instability and potential toxicity of chlorine dioxide, the preferred sterilizing gas, it is preferable 35

to transport components, which when reacted, form the chlorine dioxide gas. For example, the components may be sodium chlorite, Na₂ClO₂ and chlorine gas, Cl₂.

Appropriate control signals 7 are fed by the electronic control circuitry 100 to the valve and pump block 20 and chamber 10 for controlling components of the system. Furthermore, feedback signals 8 from the controlled components are fed back to the control circuitry 100 so that the controller can monitor the state of the system and signals 14 are coupled to display panel 110 for informing the operator of the status of the system.

Additionally, a cartridge sense signal 12 is fed from the attached gas cartridge (Cl₂ component cartridge) to indicate that a gas component cartridge has been coupled into the system.

General Functions

15

FIG. 2 shows the arrangement of valve and pump block 20 in more detail. Valve and pump block 20 includes a series of valves V1, V2, V3, V4, V4a, V5, V6, V7, V8, V9 20 and V10, pumps P1 and P2, air filter 13, a detoxifier 22 for detoxifying the evacuated chlorine dioxide gas, which may be implemented as explained in the above copending patent applications, and appropriate sources of water vapor, nitrogen, Cl, gas, air, and sodium 25 chlorite. As shown in FIG. 2, some of the valves are merely sequenced, while others are controlled in response to selected ones of the values of the measured process variables, e.g., gas concentration, humidity level and pressure. For safety reasons, each valve (V) 30 is fitted with two limit switches (LS) to indicate the open (e.g. LS2o) or closed condition (e.g. LS2c) of the valve. In the attached software listing, the open limit switches are referred to by the designation LSOx and the closed limit switches by the designation LSCx.

switches must be in their proper positions at the proper times during the entire cycle in order that the cycle not be aborted. In addition, a number of lights are provided on a display panel, as shown in FIG. 5, which indicate the progress of the sterilization cycle or the occurrence of possible fault conditions. A cycle can be started by the operator, after the chamber door 11 is closed, by momentarily depressing the START-CYCLE (S1) switch. See FIG. 1. Thereafter the cycle proceeds automatically according to a program stored in the microprocessor memory of the electronic controller 100. This process will be described in more detail below.

Furthermore, in order to provide redundancy, a number of manually controlled valves, e.g. valves V₉ and V₁₀, are provided in case valves V₃ and V₈ do not open. These valves can be manually operated by service personnel so that potentially toxic gases can be removed via detoxifier 22 in the event valves V₃ and V₈ fail to open when sterilizing gas is in the chamber. An auxiliary vacuum pump is also provided so that the gas can be drawn out via the manually operated valves.

Sterilization Cycle

The sterilization cycle is an interlocked sequence of events and consequent actions under microprocessor control. The steps of this sequence are detailed in the state diagram of FIG. 6 and state output matrix of FIG. 7. These steps are performed by a sequencing program, the flowchart for which is shown in FIGS. 7A and 7B and the details of which are disclosed in the program listing contained in the appendix. Two types of events occur during the sequence, independent and dependent events. Some independent events are external events and include contact input signals to the controller from the controlled valves (e.g., the limit switches), and are referred to by the symbols XOx to X3x in FIG. 3. Each

10

20

Safety Considerations

contact input signal is one bit of an eight bit word and the collection of such control input signals shall be referred to herein generally as digital inputs (DIN). Independent events also include the reception of signals 5 corresponding to measured or analog process values (AIN), such as pressure, temperature, humidity and sterilizing gas concentration. The measured value signals are associated with logical comparison operations performed by the controller. Other independent events are internally declared, and these typically result in the illumination of an indicator light on the display panel, shown in FIG. 5. The controller evaluates the dependent events, which are logical combinations of independent events, to single TRUE or FALSE results. When the dependent event becomes true, a corresponding action is performed, i.e. the control system moves to a new process state, defined by the state output matrix of FIG. 7. If the dependent event is not true, the controller holds the process state in its memory and waits for a period of 50 milliseconds before reevaluating the dependent event. In the case of a system failure, the system automatically transfers to an appropriate ABORT state immediately, as will be described in greater detail below. This process continues until the cycle has been completed or aborted.

The sterilization system is provided with a number of checks to insure correct operation of the various valves and other components. As will be described in more detail below, 30 interlock software implemented by the controller main timing program confirms the correct position of all valves every 6250 microseconds. An alarm condition is declared any time a valve is not in its commanded state. The operation of these interlocks differs from typical relay logic, or programmable

logic controllers, in that interlock checking continues after

valve actuation has taken place and can lead to different failure programming (ABORT states) at each process stage. correct status of a valve is latched into memory after actuation is confirmed, and this latched condition is checked 5 every 6250 microseconds. FIG. 8 is a block diagram of the safety-interlock components necessary to perform this checking routine. Failure to pass either an initial event-timeout condition following actuation or any subsequent status check will result in abnormal termination of the sterilization 10 cycle. A sequence of control actions for safe termination of the cycle is defined for every point in the sterilization cycle, and is initiated immediately in the event of any abnormal (ALARM) process condition. This intensive status checking according to the invention prevents deliberate 15 bypassing of the interlock switches in the field, since if a limit switch is bypassed, at some point in the system cycle the switch will be determined to be in an improper position, thus causing the system to enter an ABORT state.

As shown in FIG. 8, the safety interlock system includes a 20 Sequencing Program 120 stored in internal memory (ROM) of the electronic controller. Sequencing program 120 is identified in a listing of the program resident in memory attached hereto, as SEQ and the flowchart for this program is shown in FIGS. 7A and 7B. Also stored in memory is a Contact Status 25 Check program 122 and a series of masks 124 which are determined by the particular point in the sequence program. The Contact Status Check program is identified in the appendix as CSC and a flowchart therefor is shown in FIGS. 18 and 19. Inputs 126, which are images stored in memory of actual input 30 signals from both "open" limit switch contacts 127 (closed when a valve is open and open when a valve is closed) and "closed" contacts 129 (closed when a valve is closed and open when a valve is open) are provided, as well as inputs from other components, such as the sterilizing chamber 10 door 11. 35 A series of contact outputs 125 are also provided by the particular state of the sequencing program. The Contact

Status Check program 122 compares the contact inputs with the contact outputs 125. Whenever an input differs from the desired value, as established by the output, an alarm condition is declared if, and only if, a corresponding bit is turned on in the Mask 124. This safety feature detects any incorrect valve position immediately. A hardware implemented watchdog timer 132 is utilized to provide an extra level of safety by disabling all outputs to the valves 130 by opening electronic switches 134 when the timer times out if the microprocessor controller should fail, thereby preventing energization of any of the valves in the valve and pump block 20 in the event of a computer failure.

FIGS. 7A and 7B are a flowchart for the sequencing program The sequencing program is entered from another program, 15 called the Main Dispatching Program, which essentially checks for flags generated at appropriate time intervals and which determines when specific functions should be performed. As shown in FIG. 7A, when the sequencing program is entered, the current state of the system is retrieved from memory, as shown 20 at 180. The current state is stored in a register 210a in internal CPU RAM, as shown in FIG. 12. The organization of internal CPU RAM will be discussed in more detail in connection with FIG. 12 later. At 182, a check is made to determine if the state exceeds the maximum state number. If it does, an ABORT state, state 31, to be discussed in more detail in connection with FIG. 6, is entered at 184. Otherwise, the conditions for the next state are performed at 186 by entering the program ST, the flowchart for which is shown in FIG. 7B.

30 As shown in FIG. 7B, program ST first evaluates each dependent event to a single true or false result, as shown at 188 and 189. Each dependent event is a logical combination of a number of independent events, each of which must be specified if the dependent event is true. If the dependent event is not true, a hold flag (FØ) in a memory location in the

microprocessor internal RAM (see FIG. 12) is set at 190. Otherwise, the next state is set at 192 and a new ABORT state, if a new ABORT state is required, is set, but not entered, at 193.

- At 194, the timeout for the previous event must be disabled so that the timeout will not cause an alarm condition to be generated, which could cause an ABORT state to be reached. Timeouts are provided by program implemented timers, which monitor for the occurrence of a specified action, e.g. the
- 10 movement of a valve, within a preset time defined by the timer. If the specified action has occurred, the timeout must be disabled because the timer continues to run. In order to disable the timeout, as shown in FIG. 18, a flag in the Timer Counter Enable Register (TCEN) 207 in internal RAM (FIG. 12)
- is cleared. In this way, when the flag for the timer is set into the Timer Counter Flag Register (TCFL) 206 (FIG. 12) when the timer runs out, no alarm will be generated. If a timeout alarm is generated, a bit TMOF is set in the STATUS register, as shown in FIG. 18.
- 20 At 195, the masks are cleared, i.e., bits corresponding to the particular events which are to take place are set to a "don't care" condition, so that the change of the corresponding bits in the contact outputs do not set off an alarm condition by the contact status check program. At this point, the action
- 25 may be performed, as shown at 196. Subsequently, the timeout count for the action is loaded into the appropriate one of the timer registers 200 (FIG. 12) as will be explained in more detail later. The action timeout flag is then enabled to monitor for the timely occurrence of the current monitored
- 30 action as shown at 197. The hold flag FØ is then cleared at 198 and a return is made to the flowchart of FIG. 7A, to the point denoted SEQR.

At 200a, a test is performed to determine if an alarm or timeout condition has occurred. If an alarm or timeout has occurred, the current state is set to the current ABORT state at 201 immediately. Then, the hold flag FØ is checked at 202 to determine if it has been set. If it has, a return is then made to the background or main dispatching program from which all subroutines are entered. If flag FØ has not been set, the system remains in the sequencing program to continue to the next state and only exits once flag FØ is set.

FIG. 19 shows the contact status check program in more detail. As shown, the contact input status corresponding to the 10 contact inputs are stored in appropriate locations in the internal RAM of the system microprocessor. The memory locations are as indicated. See FIG. 12. The same is done for the contact output status bits, which specify the events to occur for a particular state. The Masks MSK0-MSK3, also stored in internal RAM, are evaluated by the contact status check program. If the contact inputs vary from the contact outputs, an alarm condition is generated by setting a bit in the status register 204, which is a location in RAM (see FIG. 12), but this is only done if the corresponding bit in the 20 Mask is turned on. If the bit is off, indicating that a change of the corresponding output is to be allowed to occur, no alarm will be generated, and the contact outputs will be written into an output buffer, to be described in more detail below, to actuate the appropriate controlled or sequenced 25 component, e.g., a valve or pump, without operating an alarm.

Additional safety features are also provided for in the system. As discussed above, manually actuable valves V9 and V10, operated by service personnel, and auxiliary pump P2 are provided in the event valves V8 and V3 and main pump P1 do not operate properly, thus providing a degree of redundancy. Furthermore, as shown in FIG. 2, safety features are provided to prevent the possibility of excessive temperatures and pressures in the sterilizing chamber 10. A thermally activated switch lla is provided in series with heater HTØ1 in the chamber to detect excessive temperature. For example,

should the heater HTØl fail to turn off, the thermostatic switch lla will sense an excessive temperature and interrupt the circuit.

Additionally, should excessive pressures develop in the chamber, a pressure relief valve 9 is provided for venting gases in chamber 10 through a second detoxifier 22a to the atmosphere.

Also provided is a check valve 15 in series with valve V4 which supplies sterilizing chlorine gas to the system. Check 10 valve 15 prevents the possibility of nitrogen gas from the nitrogen cannister pressurizing the chlorine gas cannister should valves V4 and V4a fail to close. Check valve 15 only allows chlorine gas to flow out of the chlorine gas cannister and prevents nitrogen gas from flowing into the chlorine gas 15 cannister if valves V4 and V4a fail to close.

Operator Interactions

The apparatus and sterilization cycle of the system according to the invention provide for minimal operator intervention and maximum safety. FIG. 5 shows an embodiment of a display panel 20 for the invention showing the various display lights. Certain lights are provided but not used, for expansion purposes. sterilization cycle cannot be initiated until the chamber 10 door 11 has been properly closed. The DOOR-OPEN light (LTI) will then be extinguished, as shown by LTO1 changing state 25 from a "1" state in state 1 to a "0" state in state 2 of FIG. 7, and the READY-FOR-CYCLE light (LT11) will be illuminated. See also FIG. 5. To start the cycle, the operator merely presses the START-CYCLE (S1) switch (see FIG. 1) when ready. Thereafter, no operator intervention is required until the cycle ends, with illumination of the REMOVE-LOAD light (LT17), or until an alarm condition has halted the cycle. latter eventuality, one of the alarm lights indicating the failure will be on. The operator notes which lights are on,

takes the necessary action and then presses the ABORT-RESET (S2) switch when ready to cycle the system back to a defined condition and to avoid the failure condition, if possible. For example, if the PURGE-FAIL light (LT5) is on, due to the possibility of an empty nitrogen tank, the tank should be replaced before pressing the S2 switch. Similarly for other failure modes, an attempt should be made to diagnose and remedy the failure condition before pressing switch S2. The subsequent actions to abort the cycle are then predetermined and automatic. No further operator intervention is necessary. Furthermore, redundancy has been provided in the system so that if a component fails, another component, e.g., a pump or valve, can take the place of the failed component so that the system can be brought out of its failure state.

15 Control Circuitry Design

The overall design of the electronic control circuitry 100 is shown in FIG. 3. The controller is microprocessor controlled, and preferably utilizes a type 8031, 8051 or 8751 microprocessor CPU 102 manufactured by Intel Corp., because of 20 the ability of these processors to perform Boolean arithmetic on bit addressable data. The CPU 102 includes self contained Random Access Memory (RAM) and Read Only Memory (ROM). Furthermore, the controller may include external ROM 104 and a non-volatile Shadow RAM (SRAM) 106 which may be a type X2210 25 manufactured by XICOR Inc. and which, as discussed heretofore, stores critical data after power-down. The controller also includes a clock crystal 108, input latch 113 receiving Digital INputs (DIN), an A/D Converter 114 and filter 114a for Analog INputs (AIN), an output latch 117 for Digital OUtputs 30 (DOU), and a WatchDog Timer 112 (WDT). The latter timer is arranged to disable all outputs to the valves to their denergized state upon failure of the microprocessor, as described above with reference to FIG. 8. Analog to digital converter 114 and analog filter 114a, convert the analog inputs from the measured gas concentration, temperature, 35

humidity and pressure parameters to digital data.

Central processor 102 is coupled to an address/data bus 116, which also couples RAM 106, ROM 104 and a bus tranceiver 105. An address latch 103 is enabled by a line 107 from the CPU/102, and latches addresses to a further bus 109, the Read/Write and Address Bus. Bus 109 allows the DIN Latch 113, A/D converter 114, a time stamp clock 119 and DOU latch 117 to be addressed at the appropriate times during execution of the sterilization sequence program, i.e., when CPU 102 calls for input data from the various valve limit switches, DIN latch 113 is addressed. At other times the A/D converter 114 or DOU latch 117 will be addressed.

Two decoders, a read enable decoder 120 and a write enable decoder 122 are coupled to bus 109 and allow latches 113 and 15 117 and A/D converter 114 to be either read from or written to. Appropriate read/write commands are coupled on lines 126, for controlling the decoders.

Furthermore, a data bus 124 is also provided for reading data from or to the input and output latches and A/D converter.

20 Several additional control lines are also employed, including a data bus enable 125 and RAM command lines 127. Line 125 enables bus transceiver 105 only for very short intervals and only during input/output (I/O) subroutines (e.g., subroutines WCO (Write Contact Outputs), RCI (Read Contact Inputs) and RAI (Read Analog Inputs), see appendix), when input and output operations are being performed, e.g., writing output information to DOU latch 117 for controlling the valves. In this way, data on the data bus 124 for actuating the various valves of the valve and pump block cannot be transmitted to the valves except under limited circumstances. This provides an additional degree of system safety. Furthermore, bus transceiver XCVR 105 is bi-directional and the direction of data transfer is controlled by one of the read and write

lines, as shown.

RAM command lines 127 issue signals to shadow RAM 106 so that failures can be logged permanently and other critical data can be stored in the event of a power failure.

A reset line 129 is also provided between the Write enable decoder 122 and watchdog timer 112 and an enable line 130 is provided between timer 112 and DOU latch 117. As previously indicated, timer 112 monitors CPU 102 for proper system operation. Normally, CPU 102 constantly resets the watchdog timer via line 129. In the event of a CPU malfunction, the reset signal will fail to appear in time and the timer 112 times out and removes the output enable signal on line 130. The removal of this signal disables all DOU latch 117 outputs, thus preventing valve energization in the event of a CPU failure. Accordingly, a still further degree of safety has been provided in the system described.

Since the elements of the controller are coupled to data buses 116 and 124, as shown in FIG. 3, they have been assigned memory addresses through which they can be accessed by the 20 microprocessor. FIG. 3A shows one arrangement of these addresses, for reference. As indicated above, certain of the devices, such as the SRAM 106 and DOU latch 117, are provided so that the data they contain can only be changed when bits of the microprocessor port lines are sequenced properly. This is a safety feature which prevents some microprocessor failure modes from causing undesired changes in memory contents or valve positions.

All processor and program timing is derived from the basic clock oscillator 108, which preferably has a frequency of 5.9904 MHz. FIG. 4 illustrates the relationship between the various frequencies used. As indicated in FIG. 3, provision may also be made to add a precision clock 119 to the system, which can be read by way of the data/address bus or via a

20

serial data communications line 118 to provide a clock-calendar for time-stamping the process data.

As shown in FIG. 4, basic processor timing is provided by the CPU internal crystal controlled clock 108. The clock 108 frequency is divided by twelve by CPU internal counter stages 130 and 132 to provide the CPU Address Latch Enable (ALE) signal of 499,200 HZ. The ALE signal is used to strobe address latch 103 so that addresses can be placed on bus 109 and further controls the operation of A/D converter 114.

10 Signal ALE is also coupled to further internal divider stages 134 and 136. Divider stage 134 provides a signal designated TIMER 1, which is further divided by an internal counter stage 138 into a 1200 bit/sec signal for serial data transfer, which optionally may be provided to transmit system data to remote 15 locations via serial line 118.

Counter stage 136 provides an interrupt, TIMERO. TIMERO provides a transition every 6250 usecs and allows the main timed function program, TMRO, to read all contact inputs and analog inputs and write all contact outputs every 6250 usecs. The operation of this program and other programs of the operating system will be described in more detail later.

The TIMERO interrupt is then further divided by program TMRO software counter stages 142, 144 and 146, to provide the respective program execution signals designated as TIC, SEC and MIN, which occur at period of 50 msecs, 1 sec and 1 min, respectively. These will be discussed in further detail below.

FIG. 3A details the assignment of addresses on address bus 109. As shown, the bus 109 is a 16 bit bus. Internal CPU RAM 30 is assigned address space 00-FF and bits A0 to A7 on the bus 109 identify the RAM locations. Internal ROM is identified by bits A_0 to A_{15} , with bits A_{12} - A_{15} always being O's, as

shown. Addresses from 0000 to 0FFF are used. The other components, external ROM 104, external RAM 106, clock 119, A/D converter 114, DIN latch 113, DOU latch 117 and watchdog timer 112 are assigned the addresses indicated in FIG. 3A. As shown, the DIN and DOU latches each are capable of latching 4 eight bit words, the DIN latches from the various limit switches and other contact inputs and the DOU latches to the various valves, pumps, etc. Digital inputs DIN and digital outputs DOU are each subdivided into four words of 8 bits each, and all eight bits of each group are accessed at one time by the respective addresses indicated in FIG. 3A.

As shown in FIG. 3, the analog pressure, temperature, humidity and chlorine dioxide gas concentration parameters are fed from respective sensors 114c to respective amplifiers 114d, e, f 15 and q. In order to provide an additional degree of system safety when sterilizing chlorine dioxide gas is being evacuated from the sterilizing chamber, it is important that the chlorine dioxide gas concentration levels be accurately measured. Accordingly, amplifier 114g for the gas 20 concentration signal is switched into a high gain mode by a control signal Y37 during the time when the sterilizing chamber is being evacuated. In this way, A/D converter 114 will compare the input concentration analog signal with a greater number of quantizing levels, thus providing a more 25 accurate indication of the actual concentration. At all other times, amplifier 114g will remain in a low gain mode. For example, when chlorine dioxide levels are being measured in the chamber for purposes of determining an adequate sterilizing concentration, much higher concentration levels 30 are being measured, and accordingly, A/D converter 114 provides an accurate digital signal corresponding to the analog concentration level. Therefore, amplifier 114g can remain in a low gain mode. Amplifier 114g may be switched to a high gain mode by signal Y37 changing from a "0" to a "1".

35 The system data-base may be thought of as being divided into

external and internal sections. The external data-base contains the Contact inputs (CCI), which are comprised of the digital inputs DIN; the Contact Outputs (CCO), which comprise the digital outputs DOU; and the Analog INputs (AIN). 5 of the external data-base are maintained in an internal data base comprising locations in internal RAM by subroutines of the TIMERO program (TMRO), which is invoked every 6250 microseconds. That is, every 6250 microseconds, all contact inputs and analog measurements are read and stored in the 10 controller internal data-base and images of the contact outputs loaded in the DOU latch. With reference to FIG. 12, which is a memory map for the internal data RAM of CPU 102, images of the contact inputs are stored as the variables CCIO through CCI3, and the filtered analog inputs are stored as the 15 variables ADIO through ADI7. The contact outputs are stored as variables CC00-CC03. Programs using the input data retrieve it only from these locations, and not from the input devices directly. Hence, the programs only operate on images of the inputs and outputs. In addition, the internal data-20 base includes a number of register banks, RBØ-RB3. number of timers 205 are provided including a 50 msec timer TICK (50 msec), a second timer TSEC (1 sec) and a minute timer TMIN (1 min). These timers provide timed function intervals for scheduling functions implemented at those intervals by the 25 system main dispatching program. The TICK timer times out after 50 msecs and sets a flag TICF in STATUS register 204 to be used by the main dispatching program to initiate all 50 msec timed functions, including a number of timers 200 in register bank RB3 which are invoked every 50 msecs, TTMx. 30 These timers are preferably invoked for monitoring timeout conditions for the system valves, for example.

The TSEC timer similarly times out after 1 sec and sets a flag (SECF) in STATUS register 204, to be used by the main dispatching program to initiate all 1 second timed functions, including a number of timers 200 in RB3 which are invoked every second, STMx. Similarly, the TMIN timer times out after

a minute and sets a flag (MINF) in STATUS register 204 to be used by the main dispatching program to initiate the 1 minute timed functions, including a number of timers 200 in RB3 which are invoked every minute, MTMx. The data memory also includes 5 registers in RB2 for keeping track of the current state and ABORT state used by the sequence program. Also included are the sequence status register 204, TCEN and TCFL registers 207 and 208, already discussed, for the timers, and a control register CTRL for enabling a control calculation to open or 10 close a valve. 4 bits of the control register, as shown, are used for controlling the four control loops of the system, corresponding to the measured temperature, humidity, pressure and gas concentration parameters. An array of bit masks 260 is provided in the internal data-base to permit "don't care" 15 conditions when comparing contact input and output status. Further descriptions of the data elements are found in the controller program source listing in the appendix to this specification.

More particularly, internal RAM of CPU 102 may be organized as The 256 (FF) memory locations are organized into 50 20 follows. msec, one sec and one minute timers in the timed function registers (memory locations 00 to 07); optional communications program registers (memory locations 08 to 0F) for controlling a receive buffer RBUF and transmit buffer TBUF; main 25 dispatching program registers (memory locations 10 to 17); timers 200 which are implemented at 50 msec, one second and one minute intervals by timers 205 (18 to 1D); (counters 1E and 1F); a status byte 204 (20); a control byte 206 (21); a timer enable byte TCEN (22); a timer flag byte TCFL (23); a 30 series of masks 260 for the inputs; (24-27); the contact output images CC00-CC03 (28-2B); contact input images CCI0 CCI3 (2C-2F); analog inputs ADIO-ADI7 (30-37); and set points for the measured process variables, such as temperature, pressure, concentration and humidity (38-3B). The remainder 35 of the internal RAM is assigned to the communications buffers (40 to 5F), the system stack (60 to 7F) and internal

microprocessor registers and storage (80 to FF), the use of which is known to those skilled in the art. Refer to Microcontroller User's Manual, published by Intel Corp., May 1982, document No. 210359-001. Although the entire system program is contained in internal ROM of the CPU 102, an external ROM may also be provided so as to allow additional programming capabilities. Alongside FIG. 12, the contents of the STATUS, CTRL, TCEN and TCFL registers by bit are shown.

State Sequence

- The progress of the sterilization cycle can be determined from the PROGRESS lights on the display panel, shown in FIG. 5.

 During a normal cycle the failure lights should never be on. Whether normal or aborted, both cycle and failure data will be maintained in a non-volatile random access memory or shadow

 15 RAM (SRAM). For example, after a designated number of cycles, e.g. three cycles, the gas cartridges will be discharged and must be replaced. The data concerning the number of cycles in which a cartridge has been used is stored in this memory. Also, after a predetermined number of cycles, or repeated

 20 failures, the system will be disabled until maintenance has been performed. This is a safety feature which cannot be bypassed in the field, and this data is also stored in the non-volatile memory.
- As discussed, FIG. 6 is a state diagram which defines the operation of the sequencing program of the sterilant system. FIG. 7 identifies the condition of the components identified in FIG. 2 as well as the display lamps shown in FIG. 5 for the various process states. The operation of the system may now be described in further detail.
- The system always begins in an initialization state, state 0.

 During this state, all output lines of the microprocessor in control circuitry 100 are set so as to initially deenergize all valves in the valve and pump block 20. After a short time

delay, valve V7 is opened to allow air into the chamber, as shown by a "1" appearing opposite VV07 for state 0 in FIG. 7. Furthermore, during this state, the control circuitry 100 stores in memory the state of all output ports of the microprocessor.

In states 0 and 1, the door to the sterilizing chamber 10 is in its open position. Once the door is closed, state 2 is entered. As indicated in FIG. 6, this means that the system is ready to begin its cycle. As further indicated in FIG. 7, in state 2, valves V1-V6 are closed, valve V7 remains open 10 and valve V8 is closed. Display lights LT1-LT6 are off, light LT11 (READY FOR CYCLE) is on and lights LT12-LT17 are off. The corresponding limit switches (LS) are in a position determined by the condition of the associated valve, e.g., for valve V2, which is closed, limit switch LS2o is open while limit switch LS2c is made. As indicated above, two limit switches are provided on each valve, one for the open position and one for the closed position, in order to insure the safety of the system. Both limit switches must be in their proper 20 position, otherwise a failure will occur.

When the door to the chamber 10 is open, the system is in state 1, once the initialization state has been passed. Accordingly, only LTl is on and the other lights are off, as shown in FIG. 7.

- Assuming the chamber door has been closed and the system is in state 2, if the START CYCLE switch S1 is pressed, the system moves to state 3. At this point valve V7 closes, as indicated by the "0" appearing in the column for state 3 in FIG. 6 and light LT12, CYCLE IN PROGRESS, turns on. As indicated in FIG.
- 30 2, valve V7 vents the chamber 10 via a filter 13 to the atmosphere when open. Thus, the flow of filtered external air into the chamber is stopped when valve V7 closes.

If the door is opened in state 2, an immediate return to state

1 is made.

Once in state 3, and, if V7 is closed, as indicated by the closed state of limit switch LS7c and open state of limit switch LS7o, state 4 will be entered. If valve V7 does not close within a certain time, as determined by a timeout implemented by one of the TIC timers TTMx in RB3 of the data memory, state 29, ABORT-1 will be entered. Furthermore, if an alarm condition occurs, such as the opening of a valve which should be closed, an alarm condition will be generated and the point of failure indicated on the display panel, indicating to the operator that a malfunction has occurred.

Once in state 3, if the chamber door is opened, the cycle will be aborted, as shown in FIG. 6.

Assuming V7 has closed and state 4 has been entered, the

chamber heater HTØl is turned on, as indicated by the "1" in

the column for state 4 opposite HTØl. If the temperature

within the chamber increases to a sufficient level within a

time-out period, state 5 can be entered. If not, ABORT-1,

state 29, is entered and a return to state 2 is thereafter

20 made when switch S2 is depressed. A safe operating

temperature is reached when temperature switch T1 (FIG. 2) is

actuated by the temperature of the atmosphere in the chamber

reaching the desired temperature. After this occurs, the

temperature in the chamber is controlled by turning the heater

on and off as required during the cycle, as indicated by the

notation "C" in the columns of FIG. 7 opposite "HTØl".

Once state 5 is entered, valve V1 is opened, in preparation for starting vacuum pump P1 so that the atmospheric contents of chamber 10 can be evacuated. Again, if valve V1 does not open within a timeout period, ABORT-1, state 29 is entered.

State 6 is entered when vacuum valve Vl opens within the timeout interval. At this point, the vacuum pump Pl is

20

started and light LT13 indicates that evacuation is in process. A timer is started which determines the length of time that the pump remains on.

Once in state 6, the chamber door 11 can no longer be opened, 5 because, at this point in the cycle, the chamber is under a vacuum.

In state 6, the pressure in the chamber is checked to determine if it has been reduced sufficiently so that it is less than or equal to a nominal value, defined as PEVAC. If the pressure is less than PEVAC, then state 7 is entered and valve VI is closed.

Should the pressure within the chamber be greater than PEVAC after the evacuation time has passed, indicating a less than adequate vacuum level, state 29 is entered. Similarly, if valve VI does not close within a specified time, state 29 is entered from state 7.

After the valve V1 has been closed in state 7, a leak-hold test is commenced in state 8. If the pressure after a leak-hold time is less than a nominal value PLEAK, state 9 is entered. If not, abort state 29 is entered.

In state 9, water vapor is allowed to enter the chamber, i.e., valve V6 is placed in a controlled open state, as indicated by "C" in FIG. 6, and a determination is made whether the humidity has reached a specified level in a certain time.

25 Should a nominal humidity HNOM not be reached within the specified time, state 30, ABORT-2, will be entered. Since evacuation has been completed, light LT13 is turned off and light LT14, which indicates a FILL IN PROGRESS, is turned on. By FILL is meant the supply of non-sterilizing gases into the chamber, e.g., steam and nitrogen gas. At this point, the system enters a new point in the state diagram wherein malfunctions allow the system to return to a different abort

state, state 30. The state of the various valves and displays for ABORT-2 (state 30) is indicated in FIG. 7.

In state 9, the humidity timer times out. If the humidity level is greater than a nominal value HNOM, state 10 is entered. Otherwise, state 30 is entered and the cycle is aborted.

In state 10, a humidity hold test is performed wherein the humidity level is monitored for a predetermined time period. If the humidity level is not maintained for the predetermined time, state 30 is entered. Otherwise, state 11 is entered. Valves V2 and V8 are opened and valve V5, along with valve V6, is then controlled on.

Valve V5 allows nitrogen to enter the system. At this point, even though valve V2 is open, chlorine dioxide cannot enter the chamber because valves V4 and V4A, which are controlled together, are closed.

In state 11, valve V2 is checked to determine that it has opened. If it has not opened withing a specified time, state 30 is entered. If valve V2 has opened in time, state 12 is entered, and valves V4 and V4A are controlled on, allowing chlorine dioxide to enter the chamber. A timer is started during which time the chlorine dioxide gas concentration levels in the chamber are measured. As explained previously, chlorine dioxide may be generated by the reaction of two components, C1₂ gas and sodium chlorite, Na₂ClO₃, on site. Chlorine gas is contained in a canister which can be coupled to the system via a connecting port, as known in the art. A container of sodium chlorite is coupled into the system between valve V2 and valve V4, as shown in FIG. 2. In state 12, LT14 is turned off and LT15, STERILIZATION IN PROGRESS, is turned on.

Once the gas concentration measured in state 12 has reached a

10

concentration greater than or equal to a nominal concentration CNOM within a preset time period, state 13 is entered. An acceptable sterilizing gas concentration might be, e.g., 1.0 mg/L to about 300 mg/L. Otherwise a new abort state, ABORT-3, 5 state 31, is entered. This new abort state is necessary because new conditions are now present in the sterilization chamber, since sterilizing chlorine dioxide gas is now present in the chamber. This requires a different set of procedures to be followed in the event of a failure, and accordingly, a new abort state is provided.

In state 13, a gas-hold test is commenced. If the gas concentration is greater than or equal to CNOM for a predetermined time period GTMR, state 14 is entered. Otherwise, state 31 is entered and the cycle is aborted.

15 In state 14, the temperature in the chamber is measured. it is greater than a minimum temperature TMIN but not higher than a maximum temperature TMAX, state 15 is entered and a sterilization timer is started. If the temperature is not adequate, state 31 is entered and an abort occurs. A typical operating temperature is approximately 30°C. 20

During state 15, sterilization is in progress. Valve V6, for humidity control, is still controlled open, and valves V4 and V4A are also controlled open. Should an alarm condition occur, e.g., if any condition changes, i.e., a valve does not 25 remain in its proper state, state 31 is entered. State 16 is entered only after a sterilization time STMR has elapsed, which typically might be several hours.

In state 16, valves V4, V4A and V6 are closed (if they do not close in the required timeout period, state 31 is entered), valve V3 is in a controlled state and valve V8 is still open. In state 17, light LT15 is turned off and light LT16 is turned on. Light LT15 is turned off when the sterilization timer has timed out and valves V4, V4A have closed. Light LT16

indicates that a purge is in progress. During state 17, the gases in the chamber are removed via valves V3 and V8 and detoxifier 22, labelled DUMP 22 in FIG. 2, which converts the chlorine dioxide into a harmless substance. The detoxification may be accomplished as explained in the above copending patent application S.N. 601,443, by passing the evacuated chlorine dioxide gas through a reducing agent, e.g., sodium thiosulfate. The detoxified gases are removed via valve V8 by vacuum pump P1. Should valves V3 and V8 fail to open within a timeout period, ABORT-3, state 31, is entered. During state 17, an evacuation timer is started which controls the amount of time during which chamber 10 is evacuated. State 18 is entered only if both valves V3 and V8 have opened in a predetermined time interval.

15 In state 18, once the evacuation timer has timed out past a time ETMR, state 19 is entered and valves V3 and V8 are closed. State 20 is entered when valves V3 and V8 close.

As shown in FIG. 6, should an alarm condition occur or should valves V3 or V8 fail to close within a specific time, state 31 20 is entered.

In state 20, valve V5 is in a controlled state. This allows nitrogen gas to enter the system as required and also prepares the system for the removal of any remaining sterilizing gases behind valve V2 via detoxifier 22 once valve V3 is reopened in state 22. In state 20, the pressure is checked. If it is greater than a maximum pressure PMAX, valve V5 is closed in state 21, turning off the nitrogen supply. If the pressure is less than PMAX, a new abort state, ABORT-4, state 32, is entered.

In state 21, valve V5 is checked to determine that it has closed within a prescribed timeout period. If it has not, state 32 is entered and the cycle is aborted. In state 22, the remaining sterilizing gases in the system are detoxified

via detoxifier 22 and reopened valves V3 and V8 and the gases removed. Once valves V3 and V8 have opened for a sufficient period of time, state 23 is entered but only if valves V3 and V8 have opened. In state 23, another timer, denoted the DESORB timer, is activated. This allows sterilizing gases which have been absorbed into the materials in the chamber to be removed or desorbed over a time period DTMR.

Should valves V3 and V8 fail to open, ABORT-5, state 33, is entered. In this circumstance, the operator will be instructed to manually activate valves V9 and/or V10 so that sterilizing gas can be removed from the system. The manually operable nature of valves V9 and V10 is indicated in FIG. 2 by a T above the valve symbols. If valves V9 and V10 are manually opened, state 33, ABORT-5 will automatically be entered.

If state 23 is successfully reached and the DESORB timer times out after a time DTMR, state 24 will be entered. At this point, valves V2, V3 and V8 are closed and a check is made to determine that these valves are closed. Then, state 25 is entered, during which a low-gas-hold test is performed. If the gas concentration is less than or equal to an acceptable value CMIN within a time period GHTM, state 26 is entered. An acceptable level of safety might be, for example, less than .5 ppm of chlorine dioxide. Otherwise, a dummy state 35 is entered, before a return is made to state 20 by operation of switch S2. This provides a delay time in which to open valves.

In state 25, the gain of amplifier 14g (See FIG. 3) is changed so that the amplifier is placed in a high gain mode during the 30 measurement of chlorine dioxide gas concentration levels during evacuation. This is indicated by the "l" in state 25 opposite GCl (gain change control). This provides more accurate measurement of concentration levels during evacuation, providing an extra degree of system safety, as

discussed previously. Also, in state 25, a counter CNT (see RB3 of FIG. 12) is decremented. This counter forces the system to cycle through states 25, 20, 21, 22, 23 and 24 via state 35 for a specified number of times determined by the 5 initial count in the counter CNTØ. Accordingly, state 35 will be entered whenever the concentration level CMIN has not been reached within time GHTM or if the counter CNT has not reached State 26 will be entered from state 25 when both the concentration is less than CMIN and CNTØ is \emptyset . 10 provided to insure system safety in the event the concentration sensor in the sterilizing chamber should fail. By going through a number of cycles via state 35, the gas concentration will be decreased, thus insuring that, even if the concentration sensor indicates the gas concentration 15 levels are below CMIN, the system will automatically cycle through a number of times necessary to reduce the concentration to acceptable safety levels. This is important, because if the concentration sensor failed and this additional safety feature was not provided, the system might indicate 20 that the gas concentration level was within acceptable levels of safety although it actually might not be.

In state 26, a counter is checked which is incremented each time the system cycles at least to step 26. If, e.g., the count is less than 3, a jump is made to state 28. If greater 25 than or equal to 3, state 27 is entered. In state 28, valve V5 is controlled on, and the count is then incremented. This allows nitrogen gas to enter the chamber.

If the cycle count is greater than or equal to 3, then state 27 is entered directly. In state 27, valves V2, V3, V4 and V8 are opened, and all remaining gas is dumped from the system and the Cl₂ gas in the cartridge is also dumped. Once sufficient time has elapsed, i.e., the Dump Hold time DHTM has elapsed, state 28 is entered. From state 28, the system enters state 37, during which the pressure in the chamber is monitored until it is within 5% of atmospheric pressure. At

this point light LT17, REMOVE LOAD, is turned on. At this point, state 38 is entered, light LT11 is turned on and actuation of switch S2 enables a return to state 1. The operator will be notified to replace the gas cartridge if the system has gone through state 27.

As indicated in FIGS. 6 and 7, after ABORT states 29 and 30 are entered, a return is made to state 2 after switch S2 is depressed and state 2 conditions are set.

In ABORT state 31 a return is made to state 20 and state 20

10 conditions are set once switch S2 is depressed. In ABORT state 32, return is made to state 19, and state 19 conditions are set. In ABORT states 33, 34, and 36, return is made to states 23, 25 and 37, respectively. If state 38 is reached, the operator receives an indication that the cycle is complete and light LT17 is turned on. To allow the chamber door to be opened, switch S2 is actuated, and state 1 is entered. If any ABORT state is reached, the appropriate failure light is illuminated. When a return is made to states 20, 23, or 26 from an ABORT state, the system then proceeds to cycle through the states which normally follow in the sequence.

General Software Functions

The sequencing program has already been described. Generally, software for the sterilization system controller is interrupt driven. Until an interrupt occurs a background task is always running via the main dispatching program. Upon interrupt, from any of several possible event sources, software control is passed to the appropriate interrupt handling program. This is illustrated in FIG. 9.

In FIG. 9, the main dispatching program 300 is shown. This program can also be found under this heading in the program listing attached hereto. Essentially, this program monitors for the occurrence of a timer flag indicating 50 msec, 1 sec

or 1 minute functions must be performed. These flags are stored in the status register (STAT) 204 of FIG. 12. When a flag occurs, the program 300 jumps to the appropriate timer program 318, 300 or 322. The timer programs are performed on a priority basis such that one minute functions are performed first and 50 msec (T50) functions last.

There are four sources of interrupting events: power-up, timer, communications, and power-down. Power-up, power-down and communications are external hardware interrupts, while the 10 timer interrupt, TMRO, is an internal hardware interrupt under program control. Except for power-up, each interrupt handling program saves the running processor context in the CPU stack before starting its task function, and the context is restored before resumption of the interrupted program. interrupt handler (TMRO) sequences all other non-interrupt programming functions. As discussed, it accomplishes this by passing one or more flags (MINF 312, SECF 314, TICF 316), signifying which of the timed tasks is to run, through the STATUS register 204 of FIG. 12. The main dispatching program 20 300 tests the flags and will cause the selected functions to be executed as shown by 318-322. This method permits further interrupt action while lower priority functions are being completed. Some of the functions performed at one minute, one second and 50 msec (TIC) intervals are as indicated in FIG. 9 25 at 318, 320 and 322, respectively. The descriptions to follow will explain the tasks to be performed under each category of interrupt event in greater detail.

Main Dispatching Program

Essentially, the main dispatching program looks for timer

30 flags and when one is found, calls the appropriate subroutine.

See FIG. 9. The main dispatching program may be found in the attached program listing.

Power Up

Upon power-up as shown at 310, the processor stack, register bank, and other functions must be initialized. This interrupt function does not require saving of the processor context.

5 Instead, previous process information will be read from the electrically reprogrammable memory SRAM 106, the clock 119 is reset and the process will resume from whichever state has been prescribed. The watchdog timer will be reset, and control will then pass back to the main dispatching program 10 300.

The power-up routine is found in the program listing under the program title INIT.

Power Fail

A power fail program is preferably implemented. One
15 embodiment for this program, as shown in FIG. 9, stores
critical memory contents at 312 into the SRAM 106, where the
data will be preserved until power is restored. The powerfail interrupt may be designed to occur whenever the 5 volt
logic line drops below 4.55 volts, and recovery to 4.75 volts
20 may be utilized for power-up. The power fail program can be
found in the attached program listing.

Communications

A communications feature (COM) may optionally be provided in the system according to the invention. The communications

25 program is activated every time a character is removed from a serial output buffer or enters a serial input buffer. The function of this program is to feed characters to the transmit buffer as they are sent and to remove characters from the receive buffer as they are received. Two FIFO queues may be provided to hold the input and output data streams. The communications program tests the input and output data streams

for the presence of termination or control characters. Flags are set in the event of termination characters. Programs, well known in the art, may be provided for processing control characters for typical serial interface devices connected to the control circuitry. For example, it may be desirable to transmit information for recording purposes over telephone lines to a printer or display device. Other programs, known in the art, can be employed to handle the standard modem control functions, e.g., RS232C commands. Hardware I/O lines may be provided for the necessary modem control signals. The communications program saves and restores the processor context.

Timed Functions

Timed functions in the controller occur on four levels as

15 follows: functions triggered by the TIMERO timer (every 6250 microseconds), functions initiated every 50 milliseconds (TICS), functions started every second, and functions which run every minute. Data is exchanged between these levels through defined data areas in the microprocessor data-base, as indicated more clearly in FIG. 10. The TMRO program also accesses the input and output devices connected to the controller. The control function (CTR), which is activated every second, transmits valve commands to the upper four bits of the CTRL register when enabled by the lower four bits of the CTRL register on a bit by bit basis, as shown in FIG. 24.

As shown in FIG. 4, timer interrupts (TMRO) occur at intervals of 6250 microseconds (6.25 milliseconds). At each interrupt, the TMRO program is entered, and all timed functions are scheduled. As the basic cycle time of the processor is approximately two microseconds, 3120 instruction cycles will elapse before the next such interrupt. Some of this time is used at each timer interrupt to perform data gathering and interlock functions, e.g. the analog inputs and data inputs are read and stored in CPU internal RAM. This is indicated at

330 in FIG. 9. Immediately following a timer interrupt the processor context will be saved in the appropriate registers. The interrupting timer, TIMERO, will then be reset and restarted. Program functions which are to occur at intervals of 50 msec., 1 sec., and 1 min. will be scheduled as shown at 332, by passing flags, as discussed, whenever the respective time interval has elapsed. Data inputs, status checks, and outputs are performed next. Finally, the previous program context is restored, and an interrupt return is executed. If any timed events are to occur, they will be performed in sequence by the main dispatching program. Otherwise the main dispatching program will be resumed.

The basic timer program, which is executed for each timer, is shown in the flowchart of FIG. 22. As shown, the timer is first decremented and a check is made to determine if the timer has timed out, i.e., reached a count of 0. If so, the corresponding timer flag is set in TCFL register 208 shown in FIG. 12. If not, the corresponding flag is cleared. Then the program is executed for the next timer, and once all timers have been completed, a return is then made to the main dispatching program.

The decrement timer function is shown in FIG. 23. As shown, when a timer is decremented, a flag is set in the TCFL register if the time has timed out, and the current count is then stored in the appropriate timer register 200.

1. TIMERO Timer (TMR0)

The lowest level timed function, occurring every 6250 microseconds, is initiated by the interrupt TIMERO. This is indicated in the uppermost portion of FIG. 10, which is a flowchart for the various timed functions. After saving the processor context, the first function of the TMRO interrupt program is to reset and restart the timer as indicated at 400. This is performed by a subroutine RRT. In FIG. 10, the

corresponding program for implementing the desired function is indicated above the flowchart symbol, and can be found in the listing in the appendix. The TMRO program is a time-critical function. Once the timer has been restarted, all of the 5 contact inputs to the controller are read into their corresponding memory images, CCIO - CCI3 as shown at 410 and These images reside in a portion of the microcomputer memory which is bit addressable. This greatly facilitates logical processing. The subroutine for implementing this 10 function is shown in FIG. 17 and is also shown in the attached program listing as subroutine RCI. The contact output information is also located in this memory, at CCOO - CCO3 and is indicated in FIG. 10 at 425. The interrupt program next performs a masked comparison of the contact input and output 15 status bits, using bit masks 415 also stored in this memory This is shown at 420. If any bits do not match their corresponding desired outputs, when masked for "don't care" conditions, an alarm condition is set by setting a bit in the STATUS register 204 (FIG. 12), as shown at 430.

Timeout alarms are also implemented by the TMRO program. A subroutine CSC2, as shown in FIG. 18 and the attached program listing, shows how timeouts are determined. When a timer times out, e.g., a timer for determining whether a valve has closed or opened in time, a flag will be set in the timer flag register TCFL. If the setting of the flag requires an abort upon failure, e.g., if the failure of a valve to close in time is to cause an abort condition, then a flag must be set in the timer enable register TCEN. This informs the timeout alarm program that an alarm condition should be set, which will cause the alarm condition to be loaded into the STATUS register. This will cause transfer to an ABORT state by the sequencing program.

Next the current contact output status is loaded from its memory image into the output contact latch by program WCO, as shown at 435. Finally, as shown at 440 and 445, the current

analog input data 445 is read (RAI), exponentially filtered (FILTER), and stored in the correct memory locations outside the bit addressable space. See FIG. 21. Eight timer interrupts take 50 milliseconds. Thus, a well-filtered analog input scan of all eight analog inputs (only four need be used for the four control loops corresponding to gas concentration, pressure, temperature and humidity) will be available each time the 50-millisecond program is entered. Therefore, every 50 msecs, the RAI program obtains 64 input samples, 8 for each channel, the eight samples for each channel then being averaged to obtain a single analog value for each channel. A return is then made to the main dispatching program. The TIMERO program is summarized in the flowchart of FIG. 11.

2. TIC Timer (T50)

· 15 The TIC functions are those which are performed every 50 milliseconds, and include the performance of the sequencing (SEQ) program. The first function performed is that of resetting the watchdog timer as shown at 500, because if this timer is not reset in time, all valve outputs will be disabled as described with reference to FIG. 8. Next, all tick timers (TTM) are decremented at 510, their counts stored at 512, and their corresponding status flags set or cleared at 520 in register TCFL 208 of FIG. 12. The setting of the timeout flags in the TCFL register 208 (See FIG. 12) also requires that the status of a corresponding bit be determined in the Timer Counter Enable Register (TCEN) 207 by the sequencing program, as shown. In this way, if the corresponding TCEN bit is not set, this informs the controller not to enter an ABORT state when the timer flag comes on. For example, when the 30 sterilization timer times out (approximately after 4 hours), an ABORT state should not be entered. For valve time-outs however, it is desired to abort if the timer times out and the valve has not opened or closed in time, and accordingly, the corresponding TCEN bit will be set by the sequencing program, 35 thus allowing an alarm to be generated. If the valve closes

in time, its corresponding TCEN bit will be disabled, and no alarm will be generated. Once the TICK timers have been decremented, the main sequencing logic 515 (SEQ), which controls the progression from one state to the next described hereinabove, is performed until it cannot progress further, due to a hold for a specified status condition not yet present. Then, the outputs are loaded into the contact output image in memory (CCO) at 530, e.g., the output data for the appropriate valves or heater to be controlled are stored in memory. Then, the TMRO program subroutine WCO will write the output images to the controlled devices on its next pass. The TIC function program is summarized in the flowchart of FIG. 13.

3. Second Timer (T1S)

15 Every second all one-second timers are decremented at step 550, the count stored at 552, and their corresponding status bits set or cleared (555). This includes the setting of flag bits (TCFL) and appropriate Timer Counter Enable bits (TCEN) depending on whether an ABORT is to occur at the occurrence of the timer flag. Finally, the control program 559 (CTR), accepting setpoints (557) from the sequencing program 515, loads the new output status for the controlled devices into the CONTROL register for subsequent loading into the contact output registers of internal RAM. During the next pass through the TMRO program, these outputs are fed to the controlled devices. As shown in FIG. 8, the timed functions occur in the order MIN, SEC and TICK. A flowchart for the one second program, TlS, is shown in FIG. 14.

As shown in FIG. 14, the first function for the one second timer program includes the clearing of the one second flag (SECF) in the STATUS register (see FIG. 12). All one second timers are then decremented, as shown in FIG. 23 and at 600 in FIG. 14. Program TIS then obtains the loop status from the sequencing program at 602, and determines if the corresponding

control bit in the CONTROL register 206 for the particular loop has been enabled at 604. Each loop corresponds to one of the four measured analog process variables, pressure, temperature, humidity and gas concentration. This is also shown in FIG. 24. As indicated, the lower four bits of the CONTROL register 206 correspond to the status of the four loops. If the loop is enabled, a value is determined by subtracting a measured input value, e.g., gas concentration or pressure, from a stored set point value from the sequencing program, as shown at 606. If this value is greater than 0, a corresponding one of the four upper bits in the CTRL register is set at 607. If the CTRL register bit is 0, then the corresponding CONTROL register bit is cleared, as shown at 608.

- 15 At 610, the program gets the next loop and repeats steps A-X for that loop. Then the next two loops are obtained and steps A-X repeated sequentially for those two loops. When all four loops have been performed, a return is made to the main dispatching program.
- 20 The interrelationship between the analog input data, set points, control register, control program (CTR), output loading program (CTL) and contact outputs CCO are shown in FIG. 24. As shown, program CTR retrieves analog input data ADI, setpoints SP and the control register (CTRL) status from memory. The new status for the control register is then determined in accordance with the flowchart of FIGS. 14 and the new status loaded into the CTRL register. Program CTL then loads the appropriate outputs for controlling the valves and heater into the appropriate contact output register in 30 memory. During the TMRO program these outputs are then coupled to the controlled devices by the program WCO. See FIGS. 10 and 20.

4. Minute Timer (TlM)

At one-minute intervals, an optional batch time clock 119 may be updated as shown at 610. This clock may be used to initiate the display of process conditions by an appropriate printing or display device. All one-minute timers are decremented at 620, and their corresponding status bits are set or cleared at 630. The TIM program is summarized in the flowchart of FIG. 15.

A sample listing of the software for the gas sterilant system according to the invention is appended below.

```
STITLE(PROGRAM FOR SC1 STERILIZATION CONTROLLER)
          医自己含含含含含含含含含含含含含含含含含含含含含含含含含含含含含含含含含
            CONSTANT DEFINITIONS
  MCHAN
            EQU
                       07H
                                             ; MAX A/D GRAN NUMBER
  CHMSK
            EQU
                       07H
                                            ; A/D CHANNEL MASK
  BNKO
            EOU
                      OOH
                                            ; RBO
  BNK1
            EQU
                       08H
                                            ;RB1
  BNK2
            EQU
                      10H
                                            ; RB2 -
  BNK3
            EOU
                      18H
                                            :RB3
  STATE
           EQU
                      26
                                          CURRENT STATE
  ABORT
          . EQU
                      R7
                                            ; ABORT STATE
  SSTA .
           EOU
                      0
                                            ;SRAM OFFSET FOR STATE
  SABO
           EQU
                      2
                                            ;SRAM OFFSET FOR ABORT
  SCNT
           EQU
                      4
                                            ;SRAM OFFSET FOR COUNT
  SMAX
           EQU
                      38
                                            ; MAX. VALID STATE
 VDLY
           EQU
                     . 8
                                            ; VALVE DELAY (400 MSEC.)
 HDLY
           EQU
                      2
                                            HEATER DELAY (2 MIN)
 TVAC
           EQU
                      30
                                           ; EVAC TIME (30 MIN)
 LKHT
           EQU
                      5
                                           ; LEAK HOLD TIME (5 MIN)
 PVAC
          EQU
                      242
                                           ; EVAC PRESSURE (95% FS)
 PRLK
                      223
                                           PRESS. LEAK LIM. (80% FS)
 HUMT
           EQU
                      30
                                           ; HUMIDIF. TIME (30 MIN.)
 HNOM
           EQU
                      207
                                           NOM. HUM. LEVEL (81% FS)
 HUMH
          EQU
                      90
                                           HUM. HOLD TIME (90 MIN.)
 TLOU
          EQU
                     Ð
                                           ; MIN. STERIL. TEMP. (0% FS); MAX. STERIL. TEMP. (100%)
 XAMT
          EQU
                     255
 CNCT
          EQU
                     15
                                           ; CONC. TIME (15 MIN.)
 CNOM
          EOU
                     64
                                           NOM. STERIL. CONC.
 CONH
          EQU
                     100
                                           GAS HOLD TIME (100 MIN)
 TSTR
          EQU
                     200
                                           ;STERIL. TIME (200 MIN)
TEVC
          EOU
                     30
                                           ; EVAC. TIME (30 MIN.)
PN2T
          EQU
                     15
                                           ;N2 PRESS. TIME (15 MIN)
DSRB
          EQU
                     30
                                           DESORB. TIME (30 MIN)
TLGH
          EQU
                     15
                                          ; LOU GAS HOLD TIME (15)
CNTM
          EQU
                     5
                                          MIN. NO. OF PURGE CYCLES
CMIN
          EQU
                     25
                                          ; MIN. CONCENTRATION (10%)
PATH
          EQU
                     12
                                          ;ATH PRESS. (5% FS)
PMAX
          EQU
                     28
                                          MAX OPER. PRESS. (11% FS)
TOMP
         EQU
                    15
                                          ; DUMP HOLD TIME (15 MIN)
PSP1
         EQU
                    60
                                          ; PRESSURE SETPOINT
TSP1
         EQU
                    60
                                          ; TEMPERATURE SETPOINT
HSP1
         EQU
                    60
                                          HUMIDITY SETPOINT
CSP1
         EQU
                    60
                                          CONCENTRATION SETPOINT
```

```
EXTERNAL DEVICE ADDRESSES
  ; *
           EXTERNAL SHADOU RAM
  SRAM
           XDATA
                      2000H
                                            ; SHADOU RAM ADDRESS
  ;
 ; *
           ANALOG INPUTS
 INO
           XDATA
                      4000H
                                            ; CHAN-O ADDRESS (PRESS.)
 IN1
           YDATA
                      6001H
                                            ; CHAN-1 ADDRESS (TEMP.)
 IN<sub>2</sub>
           XDATA
                      6002H
                                            ; CHAN-2 ADDRESS (HUM.)
 IN3
           XDATA
                      6003H
                                            ; CHAN-3 ADDRESS (CONC.)
 IN4
           XDATA
                      6004H
                                            ; CHAN-4 ADDRESS
 IN5
           XDATA
                      6005H
                                            ;CHAN-5 ADDRESS
 ING
           YDATA
                      6006H
                                            ; CHAN-6 ADDRESS
 IN7
          XDATA
                      6007H
                                           ;CHAN-7 ADDRESS
 ;
 ; *
          CLOCK PORT
 CLK
          XDATA
                     4000H
                                            CLOCK ADDRESS
 ;
 : *
          CONTACT INPUTS
XO
          XDATA
                     .OC000H
                                           ;CCI-0
                                                    ADDRESS
X1
          XDATA
                     0C001H
                                           ;CCI-1
X2
                                                    ADDRESS
          XDATA
                     0C002H
                                           ;CCI-2
X3
                                                   ADDRESS
          XDATA
                     OC003H
                                           :CCI-3
                                                   ADDRESS
;
; *
          SUITCHES
SV1
          XDATA
                     OC004H
                                           SUITCH ADDRESS
;
; *
          CONTACT OUTPUTS
YO
          XDATA
                     OECCOH
                                           :CCO-O ADDRESS
Y1
         XDATA
                     OE001H
                                           ;CCO-1 ADDRESS
Y2
         XDATA
                    OE002H
                                           ;CC0-2
Y3
                                                   ADDRESS
         XDATA
                    OE003H
                                           ;CC0-3
                                                   ADDRESS
;
; *
         UATCHDOG TIMER
UDT
         XDATA
                    CED04H
                                          : UATCHDOG RESET. ADDRESS
```

; *		A-Base all	
; * * *		******	
	DSE		
TICK	ORG	05H	:TIME COUNTERS
TSEC		1	;TICK COUNT
TMIN		1	SEC. COUNT
		1	MIN. COUNT
RPUT	ORG	OCH	SIO BUFFER POINTERS
RTAK	ds Ds	1	RCV PUT POINTER
TPUT	DS	1	RCV TAKE POINTER
TTAK	DS	1	EXAT PUT POINTER
	ORG	1	XMT TAKE POINTER
TTMO	DS	188	; IIC TIMERS
TTM1	DS	1	;TTIHER-0
	ORG	1	;TTIMER-1
STMO	DS	1AH	SECOND TIMERS
STM1	DS	1	STIMER-O
	ORG	1	;STIMER-1
MTMO	DS	1CH	; MINUTE TIMEDS
MTM1	DS	1	;MTIMER-O
	ORG	1	;MTIMER-1
CNTO	DS	1EH	COUNTERS
CNT1	DS	1	; COUNTR-O
	BSEG	•	;COUNTR-1
	ORG	20H	
STAT	DATA	20H	; INTERNAL BIT SPACE
CTRL	DATA	21H	;STATUS BYTR
TCEN	DATA	22H	CONTROL BYTE
TCFL	DATA	23H	TIMER/COUNTED ENABLES
MSKO	DATA	24H	**************************************
MSK1	DATA	25H	JULIEUT MASK DECICEDS
MSK2	DATA	26H	OULPUT HASY DECICED
MSK3	DATA	27H	OULFUL MASK DECICEDS
	ORG	28H	OUTFUL MASK PECICIES
CCOB -	DATA	28H	12MGLU I/O RITE
CC01	DATA	29H	CUTPUT PORT 1
CC02	DATA	2 A H	COTPUT PORT 1
CC03	DATA	2BH	COTPUT PORT 2
CCIO	DATA	2CH	OUTPUT PORT 3
CCI1	DATA	2DH	; INPUT PORT O
CCI2	DATA	2EH	; INPUT PORT 1
C13	DATA	2FH	INPUT PORT 2
	DSEG		; INPUT PORT 3
	ORG	30H	.1831.00
DIO	DS	1	; ANALOG DATA IMAGE
DI1	DS	1	PRESS. INPUT
DIZ	DS	1	TEMP. INPUT
DI3	DS	1	HUM. INPUT
DI4	DS	1	CONC. INPUT
DIS	DS	1	CHANNEL 4 INPUT
DI 6 DI 7	DS	1	CHANNEL 5 INPUT
	DS	1	

STP3 DS 1 ;HUM. SETPOINT CONC. SETPOINT BATCH TIME CLOCK	STPO STP1 STP2	org DS DS DS	38H 1 1	; INTERNAL DATA AREA ; PRESS. SETPOINT ; TEMP. SETPOINT
		ORG	1 3CH 1	CONC. SETPOINT

```
DATA DEFINITIONS
            [STATUS & CONTROL]
            STATUS
  TICF
            BIT
                      STAT.0
                                 TICK FLAG
  SECF
            BIT
                      STAT.1
                                  :SECOND FLAG
  MINE
            BIT
                      STAT.3
                                  ; MINUTE FLAG
  RCVF
            BIT
                      STAT.4
                                  ;RCV FLAG
  XMIF
            BIT
                     STAT.5
                                  :XMT FLAG
  TMOF
           BIT
                     STAT.6
                                  ;TIMEOUT FLAG
  ALMF
           BIT
                     STAT.7
                                  ; ALARM FLAG
           CTRL
 CENO
           BIT
                      CTRL.O
                                 PRESS. LOOP ENABLE
 CEN1
           BIT
                      CTRL.1
                                 :TEMP. LOOP ENABLE
 CEN2
                      CTRL.2
           BIT
                                 ; HUM. LOOP ENABLE
 CEN3
           BIT
                      CTRL.3
                                 CONC. LOOP ENABLE
 CTRO
           BIT
                      CTRL.4
                                 ;PRESS. LOOP OUTPUT
 CTR1
           BIT
                      CTRL.5
                                 ; TEMP. LOOP OUTPUT
 CTR<sub>2</sub>
           BIT
                      CTRL.6
                                 ; HUM. LOOP OUTPUT
 CTR3
           BIT
                      CTRL.7
                                 ; CONC. LOOP OUTPUT
           TCEN
 TENO
          BIT
                    TCEN. 0
                                 :TTO ENABLE
 TEN1
          BIT
                    TCEN.1
                                 ;TT1 ENABLE
TEN2
          BIT
                    TCEN. 2
                                 ;STO ENABLE
TEN3
          BIT
                    TCEN.3
                                 ;ST1 ENABLE
TEN4
          BIT
                    TCEN. 4
                                 ; MTO ENABLE
TEN5
          BIT
                    TCEN.5
                                ; MT1 ENABLE
TEN6
          BIT
                   TCEN. 6
                                ; MT2 ENABLE
TEN7
          BIT
                    TCEN.7
                                : MT3 ENABLE
          TCFL
TFLO
          BIT
                    TCFL.0
                                ; ITO TIMEOUT
TFL1
          BIT
                    TCFL.1
                                ;TT1 TIMEOUT
TFL2
          BIT
                    TCFL.2
                                :STO TIMEOUT
TFL3
         BIT
                    TCFL.3
                                ;ST1 TIMEOUT
TFLA
         BIT
                    TCFL.4
                                : MTO TIMEOUT
TFL5
         BIT
                   TCFL.5
                                ;MT1 TIMEOUT
TFL6
         BIT
                   TCFL.6
                                :CTO UNDERFLOW
TFL7
         BIT
                   TCFL.7
                               CT1 UNDERFLOU
```

```
DATA DEFINITIONS
            [OUTPUT PORTS]
            OPORT_O
  LT01
            BIT
                      CC00.0
                                  : DOOR-OPEN
  LT02
            BIT
                      CC00.1
                                  ; EVAC-FAIL
  LT03
            BIT
                      CC00.2
                                  ;FILL-FAIL
  LT04
            BIT
                      CC00.3
                                  :STERIL-FAIL
  LT05
            BIT
                      CC00.4
                                  ; PURGE-FAIL
  LT06
            BIT
                      CC00.5
                                  ; LOAD-UNSTERILE
  LT07
           BIT
                     CC00.6
                                  ; SPARE
 LTOS
           BIT
                     CC00.7
                                  :TEST-FAIL
           OPORT-1
 LT11
           BIT
                     CC01.0
                                  ; READY-FOR-CYCLE
 LT12
           BIT
                     CC01.1
                                  :CYCLE-IN-PROGRESS
 LT13
           BIT
                     CC01.2
                                  : EVAC-IN-PROGRESS
 LT14
           BIT
                     CC01.3
                                 ; FILL-IN-PROGRESS
 LT15
           BIT
                     CC01.4
                                 ;STERIL-IN-PROGRESS.
 LT16
           BIT
                     CC01.5
                                 ; PURGE-IN-PROGRESS
LT17
           BIT
                     CC01.6
                                 ; REMOVE-LOAD
 LT18
           BIT
                     CC01.7
                                 ; SPARE
           OPORT_2
 VV01
          BIT
                     CC02.0
                                 ; OPEN-MAIN-VAC-VALVE
VV02
          BIT
                    CC02.1
                                 ; OPEN-MAIN-GAS-VALVE
VV03
          BIT
                    CC02.2
                                 ; OPEN-MAIN-DUMP-VALVE
VV04
          BIT
                    CC02.3
                                 ; OPEN-GAS-CTRL-VALVE
VV05
          BIT
                    CC02.4
                                 ;OPEN-N2-CTRL-VALVE
VV06
          BIT
                    CC02.5
                                 ;OPEN-STEAM-CTRL-VALVE
VV07
          BIT
                    CC02.6
                                 :OPEN-BREAK-VALVE
VVOS
          BIT
                    CC02.7
                                 OPEN-DUMP-VAC-VALVE
          OPORT_3
PP01
          BIT
                    CC03.0
                                ; TURN-P1-ON
HT01
          BIT
                    CC03.1
                                ; TURN-H1-ON
SPR1
          BIT
                    CC03.2
                                :SPARE
SPR2
          BIT
                    CC03.3
                                ;SPARE
SPR3
          BIT
                    CC03.4
                                ;SPARE
SPR4
          BIT
                    CC03.5
                                :SPARE
ADZC
         BIT
                    CC03.6
                                : A/D ZERO CALIB.
LGG1
         BIT
                   CC03.7
                                ; CONC. HIGH GAIN SUITCH
```

```
DATA DEFINITIONS
            [INPUT PORTS]
                                                                    •
            IPORT_0
 LSC1
            BIT
                       CCIO.O
                                    ; V1-CLOSED
 LSC<sub>2</sub>
            BIT.
                       CCIO.1
                                    ; V2-CLOSED
 LSC<sub>3</sub>
            BIT
                       CCIO.2
                                    : V3-CLOSED
 LSC4
            BIT
                       CCIO.3
                                    ; V4-CLOSED
 LSC5
            BIT
                       CCIO.4
                                    : V5-CLOSED
 LSC<sub>6</sub>
            BIT
                       CCIO.5
                                    : V6-CLOSED
 LSC7
            BIT
                       CCIO.6
                                    : V7-CLOSED
 LSC8
            BIT
                       CCIO.7
                                    : V8-CLOSED
            IPORT_1
 LS01
            BIT
                      CCI1.0
                                    :V1-OPEN
 LS02
           BIT
                      CCI1.1
                                    :V2-OPEN
 LS<sub>0</sub>3
           BIT
                      CCI1.2
                                    :V3-OPEN
 LS04
           BIT
                      CCI1.3
                                    : V4-OPEN
 LS05
           BIT
                      CCI1.4
                                    :V5-OPEN
 LS06
           BIT
                      CCI1.5
                                    :V6-OPEN
 LS07
           BIT
                      CCI1.6
                                    :V7-OPEN
LS08
           BIT
                      CCI1.7
                                    ; V8-OPEN
           IPORT_2
DSC1
           BIT
                      CCI2.0
                                    ; DOOR-SU-CLOSED
TSC1
           BIT
                      CCI2.1
                                   ; TEMP-SW-CLOSED
SUC1
           BIT
                      CCI2.2
                                   :MAN-SU1-CLOSED
SUC2
           BIT
                      CCI2.3
                                   ; MAN-SU2-CLOSED
SIOI
           BIT
                      CCI2.4
                                   ; SPARE
SI02
           BIT
                      CCI2.5
                                   ; SPARE
SIO3
           BIT
                      CCI2.6
                                   :SPARE
SI04
           BIT
                      CC12.7
                                   :SPARE
           IPORT_3
5105
           BIT
                      CCI3.0
                                   :SPARE
S106
           BIT
                     CC13.1
                                   :SPARE
S107
           BIT
                     CC13.2
                                   ; SPARE
SI08
          BIT
                     CCI3.3
                                   ; SPARE
SID9
          BIT
                     CCI3.4
                                   :SPARE
SIIO
          BIT
                     CCI3.5
                                   SPARE
SI11
          BIT
                     CCI3.6
                                   ;SPARE
SI12
          BIT
                     CCI3.7
                                   ; SPARE
```

```
MASK BIT DEFINITIONS
           MASK-REG-0
MVC1
           BIT
                      MSKO. O
                                 ; V1-CLOSED-MASK
MVC2
           BIT
                      MSKO.1
                                 : V2-CLOSED-MASK
MVC3
           BIT
                      MSKO. 2
                                 : V3-CLOSED-MASK
HVC4
           BIT
                      MSKO.3
                                 ; V4-CLOSED-MASK
MVC5
           BIT
                      MSKO.4
                                 ; V5-CLOSED-MASK
MVC6
           BIT
                      MSKO.5
                                 ; V6-CLOSED-MASK
MVC7
          BIT
                      MSKO.6
                                 ; V7-CLOSED-MASK
MVC8
          BIT
                      MSKO.7
                                 ; V8-CLOSED-MASK
          MASK-REG-1
MV01
          BIT
                     MSK1.0
                                 ; V1-OPEN-MASK
MV02
          BIT
                     MSK1.1
                                 ; V2-OPEN-MASK
MV03
          BIT
                     MSK1.2
                                 ; V3-OPEN-MASK
MV04
          BIT
                     MSK1.3
                                 ; V4-OPEN-MASK
MV05
          BIT
                     MSK1.4
                                 ; V5-OPEN-MASK
MV06
          BIT
                     MSK1.5
                                 : V6-OPEN-MASK
MV07
          BIT
                     MSK1.6
                                .: V7-OPEN-MASK
MV08
          BIT
                     MSK1.7
                                 : V8-OPEN-MASK
          MASK-REG-2
MDC1
          BIT
                     MSK2.0
                                 : DS-CLOSED-MASK
MTC1
          BIT
                     MSK2.1
                                 ; TS-CLOSED-MASK
MSC1
          BIT
                     MSK2.2
                                 ; SU1-CLOSED-MASK
MSC2
          BIT
                     MSK2.3
                                 ; SU2-CLOSED-MASK
```

		******	***
; *	INTERRUE	T VECTORS	***********
.,	CSEG	********	*****************
•	ORG		•
RSTV:		HOOOO	· · · · · · · · · · · · · · · · · · ·
-	LJMP	INIT	RESET VECTOR
;			
	ORG	000BH	•
TINT:	LJMP	THRO	;TIMERO VEĆTOR
•			, , , , , , , , , , , , , , , , , , , ,
	ORG	0013H	
PINT:	LJMP	PURF	: PUR FAIL VECTOR
;			THE FRIE VECTOR
	ORG	001BH	
TM1V:	RETI		.TIMENI MRADAD
:			;TIMER1 VECTOR
•	ORG	0023H	
SIOV:	LJMP	SIOHND	ACRAAA
	20111	SIUAND	SERIAL DATA VECTOR
. * * * * * *	********	****	
•	POUED ES		*****************
•	POUR PA	IL HANDLER .	
•	•••••		
PURF:		******	************
PWKP:	CLR	P1.6	STORE SRAM DATA
	SETB	P1.6	; DISABLE STORE
•	RETI		; INTERRUPT RETURN
;			, Rus ante

; ****	*****	****	• • • • • • • • • • • • • • • • • • • •
;*	TIMER	INTERRUPT HANDLER	*************
*****	******	*******	
THRO:	ORG PUSH PUSH PUSH MOV CLR ACALL ACALL ACALL ACALL ACALL ACALL POP POP POP	0030H PSU ACC DPL DPH PSU, #BNKO EA RRT RCI CSC UCO RAI EA DPH DPL ACC PSU	;SAVE PROC. STATUS ;SAVE ACCUMULATOR ;SAVE DP(L) ;SAVE DP(H) ;USE RBO ;DISABLE INTERRUPTS ;RESET AND RESTART TIMERS ;READ CONTACT INPUTS ;CONTACT STATUS CHECK ;URITE CONTACT OUTPUTS ;READ ANALOG INPUTS ;RESTORE INTERRUPTS ;RESTORE DP(H) ;RESTORE DP(L) ;RESTORE ACCUMULATOR
•	RETI	- .	RESTORE PROC. STATUS :RETURN FROM TIMERO INT.

```
THRE SUBROUTINES
 RRT:
                     TRO
                                           STOP TIMERO
          MOV
                     A,#LOU(-3120+7)
                                           ; LOAD COUNT(L)
          ADD
                     A,TLO
                                           CORRECT FOR OVERRUN
          MOV
                     TLO, A
                                           ; RELOAD COUNTER(L)
          MOV
                     A, #HIGH(-3120+7)
                                           REPEAT FOR COUNT(H)
          ADDC
                     A, THO
                                           GET CORRECTED HIGH BYTE
          MOV
                     THO, A
                                           ; LOAD COUNTER(H)
          SETB
                     TRO
                                           RESTART TIMER
CLOCK:
          DJNZ
                     TICK.CLK3
                                          ; IF 50-MSEC
          MOV
                     TICK, #8
                                          : RELOAD TICK COUNT
          SETB
                     TICF
                                             SET 50-MSEC FLAG
          DJNZ
                     TSEC, CLK2
                                             IF 1-SEC
          YOU
                     TSEC, $20
                                               RELOAD TSEC COUNT
          SETB
                     SECF
                                               SET 1-SEC FLAG
          DJNZ
                     THIN, CLK1
                                               IF 1-MIN
         MOV
                    TMIN, $60
                                                 RELOAD THIN COUNT
         SETB
                    MINE
                                                 SET 1-MIN FLAG
         SJMP
                    CLK4
CLK1:
                                                 END
         CLR
                    MINF
                                               ELSE, CLR MIN. FLAG
         SJMP
                    CLX4
CLK2:
                                                 END
         CLR
                    SECF
                                             ELSE, CLR SEC. FLAG
         SJMP
                    CLK4
CLK3:
                                               END
         CLR
                    TICF
                                          ; ELSE, CLR TIC. FLAG
CLK4:
         NOP
                                             END
         RET
                                         RETURN FROM TIMER PROG.
```

RCI:	Mov	Dam	
	Mov	DPTR, #XO	POINT CONTACT INPUTS
	MOV	RO. #CCIO	POINT DATA DISCO
CI1:	CLR	~~, 54	POINT DATA-BASE INAGE
	MOVX	P1.4	FOR R4:=4 DOUNTO 0 DO ENABLE I/O
•	SETB	A. ODPTR	GET INPUT
	MOV	P1.4	· DICIDIM D.
	INC	GRO, A	DISABLE 1/0
	INC	DPTR	STORE IT IN DATA-BASE
• • •	DJNZ	RO	• EUINI NEXT INDUA
·	RET	R4,CI1	POINT NEXT STORAGE
•	KEI		: RETURN
ċsc:			; KEIUKN
	CLR	ALMF	· CI Pi Pi A · ·
	MOV	A,CCO2	CLEAR ALARM FLAG
	XRL	A,CCI1	. GET VALVE OUTBURG
	ANL	A,MSK1	CUMPARE WITH 100 Thinks
	MOV	R2,A	PICKER OFF VII DOMIN ASSOC
	MOV	A,CCO2	17076 FARILAL DECIMA
	CPL	A ,	IGEL VALVE DIITRITE
	XRL	A.CCIO	MARE CLOSED MODELS
	ANL	A, MSKO	CONFARE WITH ICC TOWNS
•	ORL	_ A,R2	THE PERSON OF A PERSON OF THE
	JZ	CSC2	, and fray, result
CSC2:	SETB	ALMF	; if hismatch
Caca;	NOP		SET ALARM FLAG
	MOV	A, TCFL	• END
	ANL	A, TCEN	GET TIMEOUTS
•	JZ	CSC4	TEST IF ENABLED
	SETB	THOF	IIF (TMO.AND.TEN)
CSC4:	SJMP	CSC5	SET TIMEOUT FLAG
CSC5:	CLR	TMOF	; END
caca:	NOP		ELSE, CLEAR TIMEOUT FLAG
_	RET	•	,
UCO:			; RETURN
aco:	MOV	DPTR, #YO	
	MOV	RO, #CCOO	POINT CONTACT OUTPUTS
C01:	MOA	R4,#4	IACAMI DATA-BICE INCA-
COI:	MOV	A, GRO	**** K4:=4 DOUNTO O DO
	CPL	λ	OZI OUIPUT DATA
	CLR	P1.4	INVERT IT FOR OUTPUT
	MOVX	ODPTR,A	• ENABLE I/O .
•	SETB	P1.4	; LOAD OUTPUT LATER
	INC	DPTR	• DISABLE I/O
	INC	RO	FOINT NEXT OUTPUT
	DJNZ	R4,C01	, TOTAL NEXT DATA
	RET		; END
		•	; RETURN

RAI:	MOV MOV CLR MOVX SETB ACALL MOV INC INC DJNZ RET	DPTR, \$ INO RO, \$ ADIO R4, \$ 8 P1.4 A, @DPTR P1.4 FILTER @RO, A DPTR RO R4, RA1	; POINT FIRST ANALOG CHAN. ; POINT FIRST ANALOG DATA ; FOR R4:=8 DOUNTO 0 DO ; ENABLE I/O ; GET ANALOG DATA ; DISABLE I/O ; FILTER ANALOG DATA ; LOAD IT INTO DATA BASE ; POINT NEXT CHANNEL ; POINT NEXT DATA ; END ; RETURN
FILTER:	MOV MUL PUSH MOV MOV MUL MOV POP ADD XCH POP ADDC RET	B, #020H AB B ACC B, #0E0H A, GRO AB R2, B B A, B A, R2 B A, B	;LOAD FILT. CONST. CB ;B,A:=0.125*R(I) ;SAVE PROD(H) ;SAVE PROD(L) ;LOAD (1-CB) CONST. ;GET X(I-1) ;B,A:=0.875*X(I-1) ;SAVE HIGH BYTE ;LOAD PROD(L) ;ADD LOU BYTES ;GET HIGH BYTE ;LOAD PROD(H) ;A,R2 IS FILTERED DATA ;RETURN

;****	SCHEDUL	ED TIME FUNCTIONS	**************
T50:	*******	***********	*******************
150:	CLR MOV	TICF	CLEAR TICK FLAG
		PSU, #BNK2	:USE RB2
	ACALL	RUT	RESET WATCHDOG TIMER
		DTT ·	DECREMENT TICK TIMERS
	ACALL	SEQ	; PERFORM SEQUENCE LOGIC
	ACALL	· CTL	LOAD CONTROL OUTPUTS
	RET		RETURN TO DISPATCHING
;			•
T1K:	CLR	SECF	CLEAR 1-SEC FLAG
•		PSU, #BNK2	USE RB2
•	ACALL	DST	DECREMENT SECOND TIMERS
	ACALL	CTR	
	RET	CIK	; PERFORM CONTROL ACTIONS
	KE1		•
;		*****	
Tim:	CLR	HINF	;CLEAR 1-MIN FLAG
	MOV	PSU, #BNK2	;USE RB2
	ACALL	UBC	;UPDATE BATCH CLOCK
	ACALL	DMT	DECREMENT MINUTE TIMERS
	RET		

RUT:	Mov Clr Clr Movx Setb Ret	DPTR, #UDT A P1.4 #DPTR, A P1.4	; POINT UATCHDOG TIMER ; CLEAR ACCUMULAȚOR ; ENABLE I/O ; RESET UATCHDOG TIMER ; DISABLE I/O
UBC:	MOV CLR XCH INC XCH INC XCH ADDC XCH RET	RO, STIME C A, GRO A A, GRO RO A, GRO A, SRO A, GRO	;POINT TIME(L) ;CLEAR CARRY ;GET TIME(L) ;INCREMENT IT ;UPDATE TIME(L) ;POINT TIME(H) ;GET TIME(H) ;PROPAGATE CARRY ;UPDATE TIME(H)

```
*************************
            CONTROL CALCULATIONS
   CIR:
            MOV
                      RO, #STPO
                                          POINT SETPOINT .
            MOV
                      R1,#ADIO
                                         POINT DATA
            CLR
                      C
            MOV
                                          CLEAR CARRY
                      A, GRO
                                         GET PRESS: SETPOINT
            SUBB
                      A, GR1
                                          SUBTRACT MEAS. PRESS.
            JNC
                      CT2
                                          ; IF MV>SP
            SETB
                      CTRO
                                            INCREASE OUTPUT
           SJMP
                     CT3
  CT2:
           CLR
                                            END
                     CTRO
  CT3:
                                         ; ELSE, DECR. OUTPUT
           NOP
           INC
                                         : END
                     RO
                                         POINT NEXT SETPOINT
           INC
                     R1
                                         POINT NEXT MEASUREMENT
           CLR
                     C
                                         CLEAR CARRY
           MOV
                     A, GRO
                                         GET TEMP. SETPOINT
           SUBB
                     A.eR1
                                         ; SUBTRACT MEAS. TEMP.
           JNC
                     CT4
                                         ; IF MV>SP
           CLR
                     CTR1
                                           DECREASE OUTPUT
          SJMP
                    CT5
 CT4:
          SETB
                                         ; END
                    CTR1
 CT5:
                                        ; ELSE, INCR. OUTPUT
          NOP
          INC
                                        ; END
                    RO
                                        POINT NEXT SETPOINT
          INC
                    R1
                                        POINT NEXT MEASUREMENT
          CLR
                    C.
          MOV
                                        CLEAR CARRY
                    A. GRO
                                        GET HUM. SETPOINT
          SUBB
                    A, eR1
                                        ; SUBTRACT HUM. MEAS.
          JNC
                    CT6
          CLR
                                        ; IF MV>SP
                    CTR2
                                        : DECREASE OUTPUT
          SJMP
                    CT7
 CT6:
         SETB
                                          END
                    CTR2
 CT7:
                                        ; ELSE, INCREASE OUTPUT
         NOP
         INC
                                        ; END
                   RO
                                       POINT NEXT SETPOINT
         INC
                   R1
                                       POINT NEXT MEASUREMENT
         CLR
                   C
                                       CLEAR CARRY
         MOV
                   A, ero
                                       GET CONC. SETPOINT
         SUBB
                   A, GR1
                                       ; SUBTRACT CONC. MEAS.
         JNC
                   CTS
         CLR
                                       ; IF MV>SP
                   CTR3
                                       ; DECREASE OUTPUT
         SJMP
                   CT9
CT8:
                                        END
         SETB
                   CTR3
CT9:
                                       ; ELSE, INCR. OUTPUT
         NOP
        RET
                                       : END
                                       ; RETURN
```

```
; *
         SOFTUARE TICK TIMERS (50 MSEC)
                   ***********************
                                        :POINT FIRST TICK TIMER
                   RO. #TTMO
DTT:
         YOK
                                        GET LAST COUNT
         VOM
                   A, GRD
                                        :IF COUNT<>0
         JZ
                   TT1
         DEC
                                        : DECREMENT ACC.
                                        ; UPDATE COUNT
         VOM
                   QRO, A
                  TT1
                                           IF NOT TIMEOUT
         JZ
                                           CLEAR FLAG
                                       •
                   TFLO
         CLR
                                             END
                   TT2
         SJMP
                                           ELSE. SET FLAG
TT1:
         SETB
                   TFLO
                                        ; END
TT2:
         NOP
                                        :POINT SECOND TICK TIMER
         YOM
                   RO, STTM1
                                        :GET LAST COUNT
         VOM
                   A. GRO
                                        :IF COUNT<>0
                    TT4
         JZ
                                          DECREMENT ACC.
         DEC
                    λ
                                           UPDATE COUNT
         MOV
                    GRO.A
                                           IF NOT TIMEOUT
                    TT4
         JZ
                                             CLEAR FLAG
         CLR
                    TFL1
                                             END
         SJMP
                    TT5
         SETB
                                           ELSE, SET FLAG
TT4:
                    TFL1
                                           END
         NOP
TT5:
                                        :RETURN
         RET
; *
         SOFTUARE SECOND TIMERS
; *
; * * * *
          ***************
                                        : POINT FIRST SEC. TIMER
         VOM
                   -RO. #STMO
 DST:
                  . A, GRO
                                        :GET LAST COUNT
         MOV
                                        :IF COUNT<>0
         JZ
                    ST1
                                           DECREMENT ACC.
         DEC
                  . A
                                           UPDATE COUNT
         VOM
                    GRO.A
                                           IF NOT TIMEOUT
          J2
                    ST1
                                             CLEAR FLAG
         CLR
                  TFL2
                                             END
          SJMP
                   ST2
                                           ELSE, SET FLAG
 ST1:
          SETB
                    TFL2
                                           END
 ST2:
       . NOP
 ÷
                                        : POINT NEXT SECOND TIMER
                    RO, #STM1
         MOV
                                        :GET LAST COUNT
         VOM
                    A. GRO
          JZ
                                        :IF COUNT<>0
                    ST4
                                           DECREMENT ACC.
          DEC
                    λ
                                           UPDATE COUNT
          MOV
                    ero, A
                                           IF NOT TIMEOUT
          JZ
                    ST4
                    TFL3
                                             CLEAR FLAG
          CLR
                   ST5
                                             END
          SJMP
ST4:
                                           ELSE, SET FLAG
                    TFL3
          SETB
                                           END
         NOP
 ST5:
                                        :RETURN
         RET
```

```
SOFTUARE MINUTE TIMERS
  DMT:
            MOV
                       RO. ANTHO
                                             POINT FIRST MIN. TIMER
            MOV
                       A, ero
                                             GET LAST COUNT
            JZ
                       HT1
                                             :IF COUNT<>0
            DEC
                       A
                                                DECREMENT ACC.
            VOM
                       QRO,A
                                                UPDATE COUNT
            JZ
                      MT1
                                                IF NOT TIMEOUT
            CLR
                      TFL4
                                                  CLEAR FLAG
            SJMP
                      MT2
                                                  END
 MT1:
           SETB
                      TFL4
                                               ELSE, SET FLAG
  MT2:
           NOP
                                               END
           MOV
                      RO, #NTM1
                                            POINT SECOND MIN. TIMER
           HOV
                      A. GRO
                                            GET LAST COUNT
           JZ
                      MT4
                                            ; IF COUNT<>0
           DEC
                      λ.
                                               DECREMENT ACC.
           MOV
                      ero, A
                                               UPDATE COUNT
           JZ
                      MT4
                                               IF NOT TIMEOUT
           CLR
                      TFL5
                                                 CLEAR FLAG
                                            ;
           SJMP
                      MTS.
                                                 END
 HT4:
           SETB
                                               ELSE, SET FLAG
 MT5:
           NOP
                                               END
           RET
                                            : RETURN
 ; *
           SOFTUARE COUNTERS
 ;****
 DCTO:
          YOM
                     RO, #CNTO
                                           ; POINT FIRST COUNTER
          YOM
                     A. GRO
                                           GET LAST COUNT
          J2
                   · DC1
                                           :IF COUNT<>0
          DEC
                     À
                                              DECREMENT ACC.
          MOV
                     GRO, A
                                              UPDATE COUNT
          JZ
                     DC1
                                              IF NOT ZERO
          CLR
                     TFL6
                                                CLEAR FLAG
          SJMP
                     DC2
                                                END
DC1:
          SETB
                     TFL6
                                              ELSE, SET FLAG
DC2:
          NOP
                                              END
          RET
                                           :RETURN
DCT1:
          YOK
                     RO, CNT1
                                           : POINT SECOND COUNTER
          YOM
                     A, GRO
                                           GET LAST COUNT
          JZ
                     DC3
                                           ; IF COUNT<>0
          DEC
                    λ
                                           ; DECREMENT ACC.
          HOV
                    ero, A
                                              UPDATE COUNT
          JZ
                    DC3
                                              IF NOT ZERO
         CLR
                    TFL7
                                                CLEAR FLAG
         SJMP
                    DC4
                                                END
DC3:
         SETB
                    TFL7
                                             ELSE, SET FLAG
DC4:
         NOP
                                             END
         RET
                                         : RETURN
```

*	CONTRO	L OUTPUTS	************
;***** CTL:	MOV ANL MOV	C,CTRO C,CENO VV05,C	GET OUTPUT-O ;ALLOW IF ENABLED ;OUTPUT TO V5
;	MOV ANL MOV	C,CTR1 C,CEN1 ET01,C	GET OUTPUT-1; ALLOW IF ENABLED; OUTPUT TO H1
;	MOV MOV	C,CTR2 C,CEN2 VV06,C	GET OUTPUT-2; ALLOW IF ENABLED . COUTPUT TO V6
	MOV ANL MOV	C,CTR3 C,CEN3 VVO4,C	GET OUTPUT-3; ALLOW IF ENABLED; OUTPUT TO V4
•	RET		

```
POUER-ON INITIALIZATION
 INIT:
           VOK
                      SP. #060H
                                            :INITIALIZE STACK POINTER
           MOV
                      PSU, 4BNKO
                                            :USE RBO
           CLR
                                            CLEAR ACCUMULATOR
           MOV
                      RO.#2
                                            POINT LOVEST RAM LOC.
          MOV
                      R1,#126
                                            :FOR R1:=126 DOUNTO 0 DO
 ILP:
          MOV
                      GRO', A
                                               CLEAR MEMORY LOC.
           INC
                     RO
                                               POINT NEXT LOCATION
           DJN2
                      R1.ILP .
                                               END
          MOV
                      TICK.#8
                                            ; INITIALIZE TICK COUNT
          MOV
                      TSEC, $20
                                            ; INITIALIZE SEC. COUNT
          MOV
                     THIN. #60
                                            ; INITIALIZE MIN. COUNT
          MOV
                     PSU, #BNK1
                                            :USE RB1
          MOV
                     RPUT, $40H
                                            ; INITIALIZE RPUT POINTER
          MOV
                     RTAK, #40H
                                            ; INITIALIZE RTAK POINTER
          MOV
                     TPUT, #50H
                                            ; INITIALIZE TPUT POINTER
          MOV
                     TTAK, $50H
                                            :INITIALIZE TTAK POINTER
          MOV
                     PSU, #BNK2
                                           :USE RB2
          HOV
                     STATE. #0
                                           ;STATE:=0
          MOV
                     ABORT, #0
                                         . ; ABORT: = 0
                     SCON, #052H
          VON
                                           ;SET SERIAL PORT BITS
          MOV
                     THOD, #061H
                                           :SET TIMER MODES
          HOV :
                     87H.400H
                                           :SET SMOD:=0 IN PCON
          YOM
                     IP. #002H
                                           SET INTERRUPT PRIORITIES
          HOV
                     IE. #096H
                                           : ENABLE INTERRUPTS
          MOV
                     TLO, #LOU(-3120)
                                           :LOAD COUNT(L)
          MOV
                     THO. #HIGH(-3120)
                                           :LOAD COUNT(H)
          MOV
                     TH1,#-13
                                           ;SET BAUD RATE (1200)
          YOM
                     A. #OFFH
                                           ;SET ACCUM. ALL 1'S
          CLR
                     P1.4
                                           ; ENABLE I/O
                     DPTR, $YO
          HOV
                                           : POINT TO OUTPUTS
          KVOM
                     edptr.a
                                          :CLEAR TO
          VOM
                     DPTR, #Y1
                                           : POINT Y1 OUTPUTS
          XVOM
                     edptr.A
                                           ; CLEAR Y1
          YOM
                     DPTR.#Y2
                                           : POINT Y2 OUTPUTS
          KVON
                     ODPTR, A
                                           :CLEAR Y2
         VOM
                     DPTR, #Y3
                                           ; POINT Y3 OUTPUTS
                     eDPTR.A
          KVOK
                                           :CLEAR Y3
          SETB
                     P1.4
                                           :DISABLE 1/0
          ACALL
                     RUT
                                           RESET WATCHDOG TIMER
         MOV
                     TIME. #0
                                           :CLEAR TIME(L)
         NOV
                     TIME+1, #0
                                           ; CLEAR TIME(H)
         SETB
                     TR1
                                           :START BAUD CLOCK
         SETB
                    TRO
                                           :START TIMER
         SJMP
                    MAIN
                                           :START MAIN PROGRAM
TEST:
         RET
                                           TEST COMPUTER FUNCTIONS
```

```
SEQUENCING PROGRAM
                                          ; REPEAT
          YOK
                     PSU, #BNK2
                                              USE RB2
          YOR
                     A,STATE
                                              GET CURRENT STATE
          ADD
                     A, #NOT(SMAX)
                                              COMPARE, MAX. STATE
          JNC
                     SQ1
                                              IF INVALID STATE
          MOV
                     A. #31
                                                TAKE STATE #31
          MOV
                     STATE, A
                                               SET STATE TO #31
                                          3
          SJMP
                     SQ2
                                                END
SQ1:
          VOM
                     A.STATE
                                             ELSE, USE CURRENT STATE
SQ2:
          NOP
                                                END
          RL .
                                             MAKE IT 4-BYTE-
          RL
                                             ADDRESS OFFSET
          YOM
                    DPTR, #JMPTBL
                                             OFFSET IN JUMP TABLE
          JMP
                    GA+DPTR
                                             PERFORM STATE
SEOR:
         YOM
                    C, ALMF
                                             GET ALARM FLAG
          ORL
                    C, THOF
                                             OR WITH TIMEOUT FLAG
          JNC
                    SQ3
                                             IF (ALM.OR.TMO)
         MOV
                    A, ABORT
                                               GET ABORT STATE
         MOV
                    STATE, A
                                               SET STATE TO ABORT
         CLR
                    FO
                                               CLEAR HOLD FLAG
SQ3:
         NOP
                                               END
         JNB
                   . FO, SEQ
                                          :UNTIL HOLD
         RET
                                          :RETURN
```

END

```
MAIN DISPATCHING PROGRAM
MAIN:
          NOP
                                           :DO FOREVER
          JNB
                     MINF, MN1
                                              IF 1-HIN TIME
          LCALL
                     TIM
                                                DO 1-MIN FUNCTIONS
MN1:
          JNB
                     SECF, MN2
                                              IF 1-SEC TIME
          LCALL
                     TIK
                                                DO 1-SEC FUNCTIONS
MN2:
          JNB
                     TICF, MN3
                                              IF TICK TIME
          LCALL
                     T50
                                                DO TICK FUNCTIONS
MN3:
          JNB
                     RCVF, MN4
                                              IF RCV TIME
          LCALL
                    RCV
                                                DO RCV FUNCTIONS
MN4:
          JNB
                     XMTF.MN5
                                             IF XMT TIME
          LCALL
                     TMX
                                                DO XMT FUNCTIONS
MN5:
          LCALL
                     TEST
                                          ; ELSE, PERFORM TESTS
          SJMP
                     MAIN
                                          ; END
GTCT:
         MOV
                     λ, $1
                                          ; READ SRAM
         RET
RCV:
         CLR
                    RCVF
                                          ; RESET RCV · FLAG ·
         RET
: TMX
         CLR
                    XMTF
                                          RESET AMT FLAG
         RET
SIOHND:
         RET
                                          ;SERIAL I/O HANDLER
$INCLUDE(STATES.SRC)
```

LJMP

STATE25

```
JMPTBL:
          LJMP
         DB
                     0
         LJMP
                     STATE1
         DB
         LJMP
                    STATE2
         DB
         LJMP
                    STATE3
         DB
         LJMP
                    STATE4
         DB :
         LJMP
                    STATE5
         DB
                    0
         LJMP
                    STATE6
         DB
         LJMP
                    STATE7
         DB
                    0
         LJMP
                    STATE8
         DB
         LJMP
                    STATE9
         DB
                    0
         LJMP
                    STATE10
         DB
         LJMP
                    STATE11
         DB
                    0
         LJMP
                    STATE12
         DB
                    0
         LJMP
                    STATE13
         DB
                    0
         LJMP
                    STATE14
         DB
         LJMP
                    STATE15
         DB
                    0
        LJMP
                    STATE16
        DB
        LJMP
                    STATE17
        DB
        LJMP
                    STATE18
        DB
        LJMP
                    STATE19
        DB
                    0
                    STATE20
        LJMP
        DB
                    0
        LJMP
                    STATE21
        DB
        LJMP
                    STATE22
        DB
        LJMP
                   STATE23
        DB
        LJMP
                   STATE24
        DB
```

DB LJMP DB	O STATE2 O STATE2 O STATE2 O STATE3
LJMP	STATE35
-	_
LJMP	STATE36
DB	0
LJMP	STATE37
DB	0
LJMP	STATE38
DB	0

```
STATEO:
            YOU
                      STATE. #1
                                            :STATE:=1
            MOV
                      ABORT, #1
                                           .; ABORT:=1
           MOV
                      STAT, #00H
                                            RESET STATUS
           MOV
                      CTRL, #00H
                                            RESET CONTROLS
           HOV
                      TCEN, #00H
                                            RESET ALARMS
           YOM
                      TCFL, #00H
                                            RESET TIMEOUT FLAGS
           MOV
                      MSKO, #00H
                                            RESET CLOSED MASKS
           MOV
                      MSK1.#00H
                                           RESET OPEN MASKS
           VON
                      MSK2, #00H
                                           RESET MISC. MASKS
           MOV
                      MSK3, #00H
                                           RESET MISC. MASKS
           MOV
                      CC00, #00H
                                           RESET ALARM LIGHTS
           MOV
                      CC01, #00H
                                           RESET RUN LICHTS
           MOV
                      CC02, #40H
                                           RESET ALL VALVES
           MOV
                      CC03, #00H
                                          RESET MISC. OUTPUTS
           CLR
                      FO
                                           CLEAR HOLD FLAG
           LJMP
                      SEOR
                                           :RETURN
 STATE1:
           JNB
                     DSC1.511
                                           ; IF DOOR CLOSED
           YOM
                      STATE, #2
                                              STATE:=2
           MOV
                      ABORT,#29
                                              ABORT:=29
           CLR
                     LT01
                                              DOOR-OPEN(OFF)
           SETB
                     LT11
                                              READY-FOR-CYCLE(ON)
           CLR
                     FO
                                              CLEAR HOLD FLAG
           SJMP
                     512
                                              END
 S11:
          SETB
                     LT01
                                           ; ELSE, DOOR-OPEN(ON)
          CLR
                     LT11
                                              READY-FOR-CYCLE(OFF)
          SETB
                     FO
                                              SET HOLD FLAG
 S12:
          NOP
                                              END
          LJMP
                     SEQR
                                          RETURN
STATE2:
          JNB
                     SUC1,521
                                          ; IF START-CYCLE(PUSHED)
          MOV
                     STATE, #3
                                          ; STATE:=3
          MOV
                     ABORT, #29
                                            ABORT:=29
          CLR
                     LT11
                                            READY-FOR-CYCLE(OFF)
          SETB
                     LT12
                                            CYCLE-IN-PROGRESS(ON)
          YOM
                     CNTO, #CNTM
                                            LOAD MIN. COUNT
          CLR
                     TFL6
                                             CLEAR COUNT FLAG
          CLR
                     MVC7
                                             CLEAR VC7 MASK
          CLR
                     HV07
                                             CLEAR VOT MASK
          CLR
                     VV07
                                          ; CLOSE-BREAK-VALVE
          MOV
                     TTMO, #VDLY
                                             LOAD TIMEOUT DELAY
          CLR
                     TFLO
                                             RESET TIMEOUT FLAG
          SETB
                     TENO
                                             ENABLE TIMEOUT ALARM
          CLR
                    FO
                                             CLEAR HOLD FLAG
          SJMP
                    S23
                                             END
S21:
          JB
                    DSC1,522
                                          ; ELSE, IF DOOR-OPEN
         Von
                    STATE, $1
                                               STATE:=1
         YOM
                    ABORT, #29
                                               λBORT:=29
         CLR
                    FO
                                               CLEAR HOLD FLAG
         SJMP
                    S23
                                               END
S22:
         SETB
                    FO
                                            ELSE, SET HOLD FLAG
S23:
         NOP
                                               END
         LJMP
                    SEOR
                                         : RETURN
```

		•	
STATE:	HOV HOV CLR SETB SETB SETB CLR SETB CLR SJMP JB HOV MOV SETB CLR	LSC7,S31 STATE,\$4 ABORT,\$29 TENO HVC7 HVO7 HT01 HTH0,\$HDLY TFL4 TEN4 F0 S33 DSC1,S32 A,ABORT STATE,A LT01 F0	; IF V7 CLOSED ; STATE:=4 ; ABORT:=29 ; CLEAR TIMEOUT ENABLE ; SET VC7 MASK ; SET VC7 MASK ; TURN HEATER ON ; LOAD HEATER TIMEOUT ; RESET TIMEOUT FLAG ; ENABLE TIMEOUT ALARM ; CLEAR HOLD FLAG ; END ; ELSE, IF DOOR OPEN ; GET ABORT STATE ; STATE:=ABORT-1 ; DOOR-OPEN(ON)
	SJMP	S33	CLEAR HOLD FLAG
S32:	SETB	FO	END
533:	NOP	20	
;	LJMP	SEOR	; ELSE, SET HOLD FLAG ; END ; RETURN

\$

```
STATE4:
           JNB
                      TSC1,541
                                          FIF HEATER ON
           NOV
                     STATE, 45
                                        ; STATE:=5
           MOV
                     ABORT, 429
                                          ; ABORT:=29
           CLR
                     TEN4
                                          : CLEAR TIMEOUT ENABLE
           SETB
                     MTC1
                                           SET TEMP SU MASK
           MOV
                     STP1.#TSP1
                                          ; LOAD TEMP. SETPOINT
           SETB
                     CEN1
                                         ; ENABLE TEMP. CONTROL
           CLR
                     MVC1
                                          : CLEAR VC1 MASK
           CLR
                     MV01
                                           CLEAR VÓ1 MASK
           SETR
                     VV01
                                          ; OPEN V1
           MOV
                     TIMO, AVDLY
                                           LOAD TIMEOUT DELAY
           CLR
                     TFLO
                                           RESET TIMEOUT FLAG
           SETB
                     TENO
                                            ENABLE TIMEOUT ALARM
           CLR
                     FO
                                          : CLEAR HOLD FLAG
           SJMP
                     S43
                                            END
                                          :
 S41:
           JB
                     DSC1,542
                                        ;ELSE, IF DOOR OPEN
          MOV
                     A. ABORT
          nov .
                                              GET ABORT STATE
                                         ;
                    STATE, A
                                         ;
                                              STATE: = ABORT-1
          SETB
                    LT01
                                              DOOR-OPEN(ON)
          CLR
                    FO
                                              CLEAR HOLD FLAG
          SJMP
                    543
                                              END
 S42:
          SETB
                    FO .
                                            ELSE, SET HOLD FLAG
 543:
          NOP
                                              END
          LJMP
                    SEOR
                                         RETURN
 STATES:
          JNB
                    LS01.S51
                                         ; IF VAC VALVE OPEN
        · MOV
                    STATE, #6
                                         ; STATE:=6
          MOV
                    ABORT, #29
                                            ABORT:=29
                                         •
          CLR
                    TENO
                                            CLEAR TIMEOUT ENABLE
                                         ;
          SETB
                    HVC1
                                           SET VC1 MASK
                                         ;
          SETB
                    MV01
                                        ; SET VO1 MASK
          SETB
                    PP01
                                           TURN P1 ON
          SETB
                    LT13
                                        ; EVAC-IN-PROGRESS(ON)
         MOV
                    MTMO, #TVAC
                                        ; LOAD EVAC TIME
         CLR
                    TFL4
                                        RESET TIMEOUT FLAG
         CLR
                    FO
                                           CLEAR HOLD FLAG
                                        :
         SJMP
                    S53
                                           END
                                        ÷
S51:
         JB ·
                    DSC1,S52
                                       ;ELSE, IF DOOR OPEN
         MOV
                   A. ABORT
                                             GET ABORT STATE
                                        :
         MOV
                   STATE, A
                                        ;
                                             STATE: = ABORT-1
         SETB
                   LT01
                                            DOOR-OPEN(ON)
                                        ÷
         CLR
                   FO
                                             CLEAR HOLD FLAG
                                        ;
         SJMP
                   S53
                                       :
                                            END
S52:
         SETB
                   FO
                                       ; ELSE, SET HOLD FLAG
         NOP
                                        .
                                           END
S53:
         LJMP
                   SEOR
                                        RETURN
STATE6:
         JNB
                   TFL4,562
                                        ; IF EVAC TIME
         CLR
                   C
                                           CLEAR CARRY
                                       ;
        YOM
                   A,ADIO
                                          GET PRESSURE
         SUBB
                   A, #PVAC
                                       ; SUBTRACT PRESS. LIMIT
         JC
                   S61
                                       : IF P.LE.PVAC
        YON
                   STATE, $7
                                          STATE:=7
                                       ï
        MOV
                   ABORT, #29
                                            ABORT:=29
                                       ;
        CLR
                   MVC1
                                            CLEAR VC1 MASK
```

S61:	CLR CLR MOV CLR SETB CLR SJMP SETB MOV MOV CLR SJMP	HVO1 VVO1 TTHO, #VDLY TFLO TENO FO S63 LTO2 A, ABORT STATE, A FO S63	:.	CLEAR VOI MASK CLOSE VI LOAD TIMEOUT RESET TIMEOUT FLAG ENABLE TIMEOUT ALARM CLEAR HOLD FLAG END END ELSE, EVAC-FAIL(ON) GET ABORT STATE STATE:=ABORT-1 CLEAR HOLD FLAG
S62:	SETB	FO		: END
S63:	NOP	FU		ELSE, SET HOLD FLAG
	· · · - -			; END
_	LJMP	SEQR		:RETURN
Ī				•

÷

```
STATE7:
             JNB
                        LSC1,571
                                             ; LF V1 CLOSED
            MOV
                       STATE, #8
                                             ; STATE:=8
            MOV
                       ABORT, #29
                                             . ABORT:=29
             CLR
                       TENO
                                                DISABLE TIMEOUT
            SETB
                       MVC1
                                             ; SET VC1 MASK
            SETB
                       MV01
                                             : SET VOI MASK
            NOV
                       MINO. &LKHI
                                             : LOAD LEAK HOLD TIME
            CLR
                       TFL4
                                             : RESET TIMEOUT FLAG
            CLR
                       FO
                                             : CLEAR HOLD FLAG
            SJMP
                       572
  S71:
                                                END
            SETB
                       FO
                                             ; ELSE, SET HOLD FLAG
            NOP
                                             : END
  S72:
            LJMP
                       SEOR
                                             RETURN
  STATES:
            JNB
                      TFL4,582
                                            ; IF LEAK HOLD TIME
            CLR
                                            ; CLEAR CARRY
            MOV
                      A,ADIO
                                               GET PRESSURE
            SUBB
                      A. #PRLK
                                               SUBTRACT LEAK LIMIT
            JC
                      581
                                               IF P.LE.PRLK
           MOV
                      STATE, #9
                                                 STATE:=9
           MOV
                      ABORT,#30
                                                 ABORT:=30
           YOM
                      STP2,#HSP1
                                              GET HUM. SETPOINT
           CLR
                      TV06
                                                 CLEAR VO6 MASK
           CLR
                      MVC6
                                                 CLEAR VC6 MASK
           SETB
                      CEN2
                                                 ENABLE HUM. LOOP (V6)
           MOV
                      NTHO, #HUMT
                                                 LOAD HUM. TIMER
           CLR
                      TFL4
                                                 RESET TIMEOUT FLAG
           CLR
                      LT13
                                                 EVAC-IN-PROGRESS(OFF)
           SETB
                      LT14
                                                 FILL-IN-PROGRESS(ON)
           CLR
                      FO
                                                 CLEAR HOLD FLAG
           SJMP
                     S83
                                                 END
 S81:
           SETB
                     LT02
                                           ; ELSE, EVAC-FAIL(ON)
          MOV
                     A.ABORT
                                                GET ABORT STATE
                                           ÷
          MOV
                     STATE, A
                                                STATE: = ABORT-1
                                           ;
          CLR
                     FO
                                           :
                                                CLEAR HOLD FLAG
          SJMP
                     583
 S82:
                                                END
          SETB
                     FO
                                           ; ELSE, SET HOLD FLAG
S83:
          NOP
                                              END
          LJMP
                     SEOR
                                           :RETURN
STATE9:
          JNB
                     TFL4,592
                                           : IF HUM. TIME
          CLR
                                           ; CLEAR CARRY
          MOV
                     A.ADI3
                                              GET HUMIDITY
          SUBB
                     A, #HNOM
                                              SUBTRACT HUM. LEVEL
          JC
                     591
                                             IF HUM.GE. HNOM
                                          ;
          MOV
                    STATE, $10
                                          ï
                                               STATE:=10
         MOV
                    ABORT, #30
                                          ;
                                               ABORT:=30
         HOV
                    MINO, & HUME
                                              LOAD HUM. HOLD TIMER
         CLR
                    TFL4
                                               RESET TIMEOUT FLAG
         CLR
                    FO
                                               CLEAR HOLD FLAG
         SJMP
                    583
591:
                                               END
         SETB
                    LT03
                                          : ELSE, FILL-FAIL(ON)
: GET ABORT STATE
         YOM
                    A, ABORT
         MOV
                    STATE, A
                                               STATE: = ABORT-2
                                          ;
         SJMP
                    S93
```

END

592:	CETO	50	
593:	SETB	£0	; else, set hold flag
373:	NOP		; END
	LJMP	SEOR	: RETURN
STATE1	O: JNB	TFL4,5101	; IF HUM. HOLD TIME
•	MOV	STATE, \$11	• STATE-#11
	MOV	ABORT, #30	; ABORT:=30
	CLR	HVC2	CLEAR VC2 MASK
•	CLR	MV02	
* .	SETB	VV02	; CLEAR VOZ MASK
	CLR	MVC8	; OPEN V2
	CLR	MVOS	; CLEAR VC8 MASK
	SETB		: CLEAR VOB MASK
		VV08	; OPEN V8
	HOV	TIMO, #VDLY	' ; LOAD VALVE TIMEOUT
	CLR	TFLO .	: RESET TIMEOUT FLAG
	SETB	TENO	; ENABLE TIMEOUT ALARM
	MOV	STPO, #PSP1	; GET PRESS. SETPOINT
	CLR	MV05	; CLEAR VOS MASK
	CLR	MVC5	; CLEAR VC5 MASK
•	SETB	CENO	: ENABLE PRESS. LOOP (V5)
	CLR	FO	; CLEAR HOLD'FLAG
	SJMP	S102	; END
S101:	SETB	FO	
5102:	NOP		ELSE, SET HOLD FLAG
	LJMP	SEOR	END
:			; return

```
STATE11: MOV
                       C,LSO2
                                             ;TEST V2 OPEN-
            ANL
                       C,LSO8
                                             ; AND V8 OPEN
            JNC
                       S111
                                             ; IF (V2.AND.V8) OPEN
            HOV
                       STATE, $12
                                                STATE:=12
            MOV
                       ABORT, #31
                                                ABORT:=31
            CLR
                       TENO
                                                DISABLE TIMEOUT FLAG
            SETB
                       HVC2
                                                SET VC2 MASK
            SETB
                       MV02
                                             SET VOZ MASK
            SETB
                       MVC8
            SETB
                       MV08
                                                SET VOS MASK
                                             . .
            HOV
                       STP3, #CSP1
                                                GET CONC. SETPOINT
            CLR
                       HV04
                                                CLEAR VO4 MASK
            CLR .
                       HVC4
                                                CLEAR VC4 MASK
            SETB
                       CEN3
                                                ENABLE CONC. LOOP (V4)
            HOV
                       HTHO, #CNCT
                                               · LOAD CONC. TIMER
            CLR
                       TFL4
                                                RESET TIMEOUT FLAG
            CLR
                       LT14
                                                FILL-IN-PROGRESS(OFF)
                                             ŧ
            SETB
                       LT15
                                                STERIL-IN-PROGRESS(ON)
                                             ÷
            CLR
                       FD
                                                CLEAR HOLD FLAG
            SJMP
                       5112
                                             ŧ
                                                END
  S111:
            SETB
                       FO
                                                ELSE, SET HOLD FLAG
 S112:
           NOP
                                                  END
           LJMP
                      SEOR
                                             ; RETURN
 STATE12: JNB
                      TFL4,S122
                                            ; IF CONC. TIME
           CLR
                      C
                                               CLEAR CARRY
           MOV
                      A.ADI3
                                               GET CONC.
           SUBB
                      A, #CNOM
                                                SUBTRACT CONC. LEVEL
           JC
                      S121
                                                IF CONC.GE.CNOM
           YOM
                      STATE, $13
                                                 STATE:=13
           YOM
                      ABORT. #31
                                                 ABORT:=31
                                            i
           NOV
                      MTMO, # CONH
                                                 LOAD CONC. HOLD TIMER
                                            ;
           CLR
                      TFL4
                                                 RESET TIMEOUT FLAG
           CLR
                      FO
                                                 CLEAR HOLD FLAG
           SJMP
                      5123
                                                 END
 S121:
           SETB
                      LT04
                                               ELSE, STERIL-FAIL(ON)
           MOV
                      A, ABORT
                                                 GET ABORT STATE
           MOV
                      STATE, A
                                                 STATE: = ABORT-3
           CLR
                      FO
                                                 CLEAR HOLD FLAG
           SJMP
                      S123
                                                 END
 5122:
           SETB
                      FO
                                            ELSE, SET HOLD FLAG
 S123:
           NOP
                                               END
           LJMP
                      SEQR
                                            :RETURN
 STATE13: JNB
                     TFL4,S132
                                            ; IF GAS HOLD TIME
          CLR
                                              CLEAR CARRY
          HOV
                     A, ADI3
                                              GET CONC.
          SUBB
                     A, #CNOM
                                              SUBTRACT CONC. LEVEL
                                           ï
          JC
                     S131
                                              IF CONC.GE.CNOM
                                           •
          HOV
                     STATE, $14
                                                STATE:=14
                                           ;
          MOV
                     ABORT, #31
                                                ABORT:=31
          CLR
                     FO
                                                CLEAR HOLD FLAG
          SJMP
                     5133
                                                END
S131:
          SETB
                     LT04
                                              ELSE, STERIL-FAIL (ON)
          YOR
                     A. ABORT
                                                GET ABORT STATE
         HOV
                     STATE, A
                                                STATE: = ABORT-3
         CLR
                    FO
                                                CLEAR HOLD FLAG
         SJMP
                    S133
                                                END
5132:
         SETB
                    FO
                                           ELSE, SET HOLD FLAG
5133:
         NOP
                                              END .
         LJMP
                    SEOR
                                          ; RETURN
```

```
STATE14: CLR
                       C
                                             CLEAR CARRY
            MOV
                       A.ADI1
                                             GET TEMP.
            SUBB
                       A. #TLOU
                                            SUBTRACT MIN. TEMP.
            JC
                       S141
                                            ; IF TEMP.GE.THIN
            CLR
                                            : CLEAR CARRY
            MOV
                       A, $TMAX ···
                                          . GET MAX. TEMP LEVEL
            SUBB
                       A, ADI1
                                         . ; SUBTRACT TEMP.
            JC
                      5141
                                           ; IF TEMP LE. THAY
           YOM
                      STATE, $15
                                            .i.
           Von
                                               STATE: =15
                      ABORT, #31
                                                ABORT:=31
                                     YOU
                      MIMO, #ISTR
                                            ; -
                                                LOAD STERIL. TIMER
           CLR
                      TFLA
                                               RESET TIMEOUT FLAG
                                            ;
           CLR
                      FO
                                            CLEAR HOLD FLAG
           SJMP
                      S142
 S141:
           SETB
                      LT04
                                           ;ELSE, STERIL-FAIL(ON)
           YOK
                      A. ABORT
                                            ; GET ABORT STATE
           MOV
                      STATE.A
                                               STATE:=ABORT-3
           CLR
                      FO
                                           : CLE
                                               CLEAR HOLD FLAG
 S142:
           NOP
           LJMP
                      SEOR
                                           :RETURN
 STATE15: JNB
                     TFL4,S151
                                           ; IF STERIL. TIME
           YOM
                                          ; STATE:=16
; ABORT:=31
; PRESS. LOOP (OFF)
; HUM. LOOP (OFF)
; GAS LOOP (OFF)
                     STATE, $16
          NOV
                     ABORT,#31
          CLR
                     CENO
          CLR
                     CEN2
          CLR
                     CEN3
          CLR
                     CTRO
                                           ; PRESS. OUTPUT (OFF)
; HUM. OUTPUT (OFF)
          CLR
                     CTR2
          CLR
                     CTR3
                                              GAS OUTPUT (OFF)
          CLR
                     VV06
                                             CLOSE V6
                                           4
          CLR
                     VV05
                                              CLOSE V5
          CLR
                     VV04
                                              CLOSE V4
                                          :
          VON
                     TIMO, #VDLY
                                              LOAD TIMEOUT DELAY
                                          ;
          CLR
                     TFLO
                                          RESET TIMEOUT FLAG
          SETB
                     TENO
                                          : ENABLE TIMEOUT ALARM
          CLR
                    FO
                                             CLEAR HOLD FLAG
          SJMP
                    S152
S151:
                                             END
          SETB
                    FO
                                          ; ELSE, SET HOLD FLAG
S152:
          NOP
                                          : END
         LJMP
                    SEOR
                                          : RETURN
STATE16: MOV
                    C.LSC4
                                        TEST V4 CLOSED
         ANL
                    C, LSC5
                                        ; AND V5 CLOSED
         ANL
                    C, LSC6
                                          :AND V6 CLOSED
         JNC
                    S161
                                         :IF (V4,V5,& V6) CLOSED
         VOM
                    STATE, $17
                                         ; STATE:=17
         VOM
                    ABORT,#31
                                            ABORT:=31
                                         2
         CLR
                    TENO
                                            DISABLE TIMEOUT ALARM
                                         ;
         SETB
                    HVC4
                                         ; SET VC4 MASK
; SET VO4 MASK
; SET VC5 MASK
         SETB
                   HV04
                                      .
         SETB
                   MVC5.
         SETB
                   MV05
                                         ;
                                            SET VOS MASK
        SETB
                   MVC6
                                         •
                                            SET VC6 MASK
        SETB
                   MV06
                                         :
                                            SET VO6 MASK
```

```
CLR
                       MVC3
                                             ; ! CLEAR VC3 MASK
            CLR
                       MVO3
            SETB
                                                CLEAR VOS MASK
                       VV03
            CLR
                                                OPEN V3
                       LT15
                                                STERIL-IN-PROGRESS(OFF)
           SETB
                       LT16
                                                PURGE-IN-PROGRESS(ON)
           MOV
                      TIMO, #VDLY
                                               LOAD TIMEOUT DELAY
           CLR
                      TFLO
                                               RESET TIMEOUT FLAG
           SETB
                      TENO
                                               ENABLE, TIMEOUT ALARM
           CLR
                      FO
                                               CLEAR HOLD FLAG
           SJMP
                      S162
 S161:
           SETB
                                               END
                      FO
 S162:
                                            ; ELSE, SET HOLD FLAG
           NOP
           LJMP
                                               END
                      SEOR
                                            ; RETURN
 STATE17: MOV
                      C, LS03
                                            ; TEST V3 OPEN-
          ANL
                     C,LSO8
                                           ; AND V8 OPEN
          JNC
                     S171
                                           ; IF (V3.AND.V8) OPEN
          MOV
                     STATE, $18
          MOV
                                              STATE:=18
                     ABORT,#31
          CLR
                                              ABORT:=31
                     TENO
                                              DISABLE TIMEOUT ALARM
          SETB
                     MVC3
                                              SET VC3 MASK
          SETB
                     MV03
                                              SET VO3 MASK
          SETB
                     MVC8
                                              SET VC8 MASK
          SETB
                                          . ;
                     MVOB
          MOV .
                                              SET VOS MASK
                     MTMO, #TEVC
                                              LOAD EVAC. TIMER
          CLR
                     TFL4
                                             RESET TIMEOUT FLAG
          CLR
                     FO
                                             CLEAR HOLD FLAG
         SJMP
                    S172
S171:
         SETB
                                             END
                    FO
5172:
                                          ; ELSE, SET HOLD FLAG
         NOP
                                             END
         LJMP
                                          :
                    SEOR
                                          ; RETURN
```

```
STATE18: JNB
                       TFL4,5181
                                           : IF EVAC. TIME
            MOV
                       STATE, $19
                                                STATE:=19
            MOV
                       ABORT,#31
                                                ABORT:=31
            CLR
                       MVC3
                                                CLEAR VC3 MASK
            CLR
                       MVO3
                                                CLEAR VO3 MASK
            CLR
                       VV03
                                             : CLOSE V3
: CLEAR VC8 MASK
            CLR
                       HVC8
            CLR
                       BOVE
                                                CLEAR VOS MASK
            CLR
                       VVOS
                                                CLOSE V8
            MOV
                       TIMO, #VDLY
                                           . :
                                                LOAD VALVE TIMER
            CLR .
                       TFLO -
                                             • .
                                                RESET TIMEOUT FLAG
            SETB
                       TENO
                                                ENABLE TIMEOUT ALARM
            CLR
                       FO
                                                CLEAR HOLD FLAG
                                             1
            SJMP
                       S182
                                                END
                                            . 2
  5181:
            SETB
                       FO
                                             ; ELSE, SET HOLD FLAG
  5182:
            NOP
                                             : END
            LJMP
                       SEOR
                                             :RETURN
  STATE19: MOV
                       C, LSC3
                                             :TEST V3 CLOSED-
            ANL
                       C, LSC8
                                             ; AND V8 CLOSED
            JNC
                       S191
                                             ; IF (V3.AND.V8) CLOSED
           MOV
                       STATE, #20
                                                STATE: =20
                                            5
           MOV
                       ABORT, #32
                                                ABORT:=32
                                            . :
           CLR
                      TENO
                                               DISABLE TIMEOUT ALARM
                                            :
           SETB
                      MVC8
                                               SET VC8 MASK
                                            :
           SETB
                      MV08
                                               SET VOS MASK
                                            ;
           MOV
                      STPO, PSP1
                                               GET PRESS. SETPOINT
                                            :
           CLR
                      MYOS
                                               CLEAR VOS MASK
                                            ;
           CLR
                      MVC5
                                               CLEAR VC5 MASK
                                            ;
           SETB
                      CENO
                                               ENABLE PRESS. CONTROL (VB)
                                            .
           YOM
                      MTMO, #PN2T
                                               LOAD N2 PRESS. TIMER
                                            :
           CLR
                      TFL4
                                            :
                                               RESET TIMEOUT FLAG
           CLR
                      FO
                                               CLEAR HOLD FLAG
           SJMP
                      S192
                                               END
 S191:
           SETB
                      FO
                                            ; ELSE, SET HOLD FLAG
 S192:
           NOP
                                               END
           LJMP
                      SEOR
                                            ; RETURN
 STATE20: JNB
                      TFL4,5202
                                            ; IF REPRESS. TIME
           CLR
                                               CLEAR CARRY
           MOV
                      A, #PMAX
                                               GET MIN.PRESS. LEVEL
           SUBB
                      A, ADIO
                                               SUBTRACT PRESSURE
           JC
                      S201
                                               IF PRESS.GE.PMAX
          NOV
                      STATE, $21
                                            ;
                                                 STATE:=21
          MOV
                      ABORT,#32
                                            ;
                                                 ABORT:=32
          CLR
                      CENO
                                                 N2 LOOP (OFF)
                                           :
          CLR
                     CTRO
                                              N2 OUTPUT (OFF)
          CLR
                     VV05
                                                 CLOSE N2 VALVE
          MOV
                     TIMO, AVDLY
                                               LOAD VALVE TIMEOUT.
          CLR
                     TFLO
                                                 RESET TIMEOUT FLAG
          SETB
                     TENO
                                                ENABLE TIMEOUT ALARM
          CLR
                     FO
                                                CLEAR HOLD FLAG
                                           ;
          SJMP
                     S202
                                                END
                                           •
S201:
          SETB
                     LT05
                                           ;
                                              ELSE, PURGE-FAIL (ON)
          VOM
                     A, ABORT
                                                GET ABORT STATE
         YOM
                     STATE, A
                                                STATE: = ABORT-4
          CLR
                     FO
                                                CLEAR HOLD TIMER
         SJMP
                     S203
                                                END
S202:
         SETB
                     FO .
                                           ; ELSE, SET HOLD TIMER
S203:
         NOP
                                           : END
         LJMP
                    SEQR
                                           :RETURN
;
```

```
STATE21: JNB
                      LSC5, S211
                                           ; IF V5 CLOSED
           MOV
                      STATE, #22
                                           ; STATE:=22
           MOV
                      ABORT,#33
                                              ABORT:=33
                                           ï
           CLR
                      TEND
                                              DISABLE TIMEOUT ALARM
                                           ;
           SETB
                                         MVC5
                                              SET VC5 MASK
           SETB
                      MV05
                                              SET VOS MASK
                                           3
           CLR
                     . MVC3
                                              CLEAR VC3 MASK
                                           i
           CLR
                      EOVII
                                              CLEAR VO3 MASK
                                           :
           SETB
                      VV03
                                              OPEN V3
                                           :
           CLR
                      HVC8
                                              CLEAR VC8 MASK
                                           :
           CLR
                      HV08
                                              CLEAR VOS MASK
                                           :
           SETB
                     VVOS
                                              OPEN V8
                                           i
           MOV
                     TIMO, #VDLY
                                           •
                                              LOAD TIMEOUT DELAY
           CLR
                     TFLO
                                              RESET TIMEOUT FLAG
                                           :
           SETB
                     TENO
                                           : . ENABLE TIMEOUT ALARM
           CLR
                     FO
                                              CLEAR HOLD FLAG
           SJMP
                     ·S212
                                              END
 S211:
           SETB
                     FO
                                          :ELSE, SET HOLD FLAG
 S212:
          MOP
                                           2
                                             END
          LJMP
                     SEOR
                                         .:RETURN
 STATE22: MOV
                    C, LS03
                                          :TEST V3 OPEN-
          ANL
                     C, LSO8
                                          :AND V8 OPEN
          JNC
                     S221
                                          :IF (V3.AND.V8) OPEN
          MOV
                     STATE, #23
                                          : STATE:=23
          MOV
                     ABORT.#33
                                              ABORT:=33
                                          ;
          CLR
                     TENO
                                              DISABLE TIMEOUT ALARM
                                          •
          SETB
                     MVC3
                                             SET VC3 MASK
SET VC3 MASK
SET VC8 MASK
SET VO8 MASK
                                          ;
          SETB
                    EOVM
                                          :
          SETB
                    MVC8
                                          :
          SETB
                    MV08
                                          ;
          MOV
                    MTMO, #DSRB
                                             LOAD DESORB TIMER
                                          ;
          CLR
                    TFL4
                                             RESET TIMEOUT FLAG
                                          :
          CLR
                    FO
                                             CLEAR HOLD FLAG
                                         :
          SJMP
                    S222
                                          ;
                                             END
5221:
          SETB -
                    FO
                                             ELSE, SET HOLD FLAG
                                          :
S222:
          NOP
                                               END
                                         . ;
          LJMP
                    SEOR
                                          ; RETURN
STATE23: JNB
                    TFL4,5231
                                        ; IF DESORB TIME
          VOM
                    STATE, #24
                                         ; STATE:=24
          YON
                    ABORT, #34
                                         : ABORT:=34
         CLR
                    HVC3
                                         CLEAR VC3 MASK
         CLR
                    MVO3
         CLR
                    VV03
                                         : CLOSE V3
         CLR
                    MVC8
                                         ; CLEAR VC8 MASK
         CLR
                    BOVM
                                         ; CLEAR VOS MASK
         CLR
                    VVOB
                                         : CLOSE V8
         CLR
                    HVC2
                                         ; CLEAR VC2 MASK
         CLR
                                         CLEAR VOZ MASK
                    MVO2
         CLR
                    VVD2
                                            CLOSE V2
                                         :
         YOM
                    TIMO, #VDLY
                                            LOAD TIMEOUT DELAY
                                         ;
         CLR
                    TFLO
                                            RESET TIMEOUT FLAG
                                         :
         SETB
                    TENO
                                            ENABLE TIMEOUT ALARM
         CLR
                    FO
                                         : CLEAR HOLD FLAG
         SJMP
                    S232
5231:
         SETB
                                         : END
                    FO
S232:
                                         ; ELSE, SET HOLD FLAG
         NOP
         LJMP
                                         : END
                    SEOR
÷
                                         :RETURN
```

```
STATE24: HOV
                       C, LSC3
                                             :TEST V3 CLOSED-
            ANL
                       C, LSC8
                                           ... ; AND V8 CLOSED-
                       C, LSC2
            ANL
                                             ; AND V2 CLOSED
            JNC
                       5241
                                             ; IF (V2, V3, V8 CLOSED)
            YON
                       STATE, #25
                                                STATE:=25
                                             3
           YOK
                       ABORT,#35
                                                ABORT:=35
                                             1
            CLR
                       TEND .
                                             •
                                                DISABLE TIMEOUT MASK-
           SETB
                      HVC3
                                                SET VC3 MASK
SET VO3 MASK
                                             ;
           SETB .
                      MV03
                                             i
           SETB
                      MVC8
                                                SET VC& MASK
                                             3
           SETB
                      MVO8
                                                SET VOS HASK
                                            ;
           SETB
                      MVC2
                                                SET VC2 MASK
                                             i
           SETB
                      HV02
                                                SET VOZ MASK
           SETB
                      LGG1
                                            •
                                                SUITCH TO HIGH GAIN
           MOV
                      MTMO, #TLGH
                                            : . START LOW GAS HOLD
           LCALL
                      DCTO
                                                DECREMENT PURGE COUNT
           CLR
                      TFL4
                                               RESET TIMEOUT FLAG
           CLR
                      FO
                                                CLEAR HOLD FLAG
           SJMP
                      S242
                                               END
 5241:
           SETB
                      FO
                                            ; ELSE, SET HOLD FLAG
 5242:
           NOP
                                            ; END
           LJMP
                      SEQR
                                            ; RETURN
 STATE25: JNB
                      TFL4,S252 .
                                            : IF LOW-HOLD TIME
           CLR
                     . C
                                            ; CLEAR CARRY
           MOV
                      A, #CMIN
                                            GET MAX. LEVEL SUBTRACT CONC.
           SUBB
                      A,ADI3
           ORL
                      C,/TFL6
                                            ; OR CARRY WITH COUNT FLAG
           JC
                      S251
                                           : IF (CONC.LE.CHIN).AND.TFL4=1
          MOV
                     STATE, $26
                                                STATE:=26
          MOV
                     ABORT, #36
                                                 ABORT:=36
          CLR
                     CEN1
                                                 DISABLE TEMP. CTRL
          CLR
                     PP01
                                                 TURN PUMP OFF
          CLR
                     HT01
                                                 TURN HEATER OFF
          CLR
                     LGG1
                                                 SET LOU GAIN
          CLR
                     FO
                                                 CLEAR HOLD FLAG
          SJMP
                     S253
                                                 END
$251:
          MOV
                     A, ABORT .
                                               ELSE, GET ABORT STATE
          MOV
                     STATE, A
                                                 STATE:=35
          CLR
                     FO
                                                 CLEAR HOLD FLAG
          SJMP
                     S253
                                                 END
S252:
          SETB
                     FO
                                           ; ELSE, SET HOLD FLAG
$253:
          NOP
                                              END
          LJMP
                     SEOR
                                           :RETURN
STATE26: LCALL
                     GTCT
                                           ; DECREMENT & GET CYCLE CHT
          JNZ
                     S261
                                           ; IF LAST RUN
         MOV
                     STATE, #27
                                              STATE:=27
         MOV
                     ABORT, $36
                                              ABORT:=36
         CLR
                     HVC2
                                              CLEAR VC2 MASK
         CLR
                    MVO2
                                           CLEAR VOZ MASK
         SETB
                     VV02
                                              OPEN V2
         CLR
                    MVC3
                                              CLEAR VC3 MASK
         CLR
                    EOVM
                                              CLEAR VOS MASK
         SETB
                    VV03
                                              OPEN V3
```

S261: S262:	CLR CLR SETB CLR SETB MOV CLR CLR SJMP MOV CLR NOP LJMP	MVC4 HVO4 VVO4 HVC8 HVO8 VVO8 HTMO, \$TDMP TFL4 FO S262 STATE, \$28 ABORT, \$36 FO SEQR	; CLEAR VC4 MASK ; CLEAR VO4 MASK ; OPEN V4 ; CLEAR VC8 MASK ; CLEAR VO8 MASK ; OPEN V8 ; LOAD DUMP TIMER ; RESET TIMEOUT FLAG ; CLEAR HOLD FLAG ; END ;ELSE, STATE:=28 ; ABORT:=36 ; CLEAR HOLD FLAG ; END ;RETURN
----------------	---	---	---

```
STATE27: JNB
                        TFL4,5271
                                              : IF DUMP-TIME
             MOV
                        STATE, #28
                                              : STATE:=28
             MOV
                        ABORT, #36
                                                 ABORT: #36
             CLR
                        VV02
                                                CLOSE V2
             CLR
                        VVD3
             CLR
                                              ;
                                                 CLOSE V3
                        VV04
                                                 CLOSE V4
                                              •
             CLR
                        VVOS
             CLR .
                                                 CLOSE VB
                        FO
                                                 CLEAR HOLD FLAG
             SJMP
                        S272
   S271:
                                                 END
             SETB
                        FO
                                              ; ELSE, SET HOLD FLAG
   S272:
             NOP
                                                 END
             LJMP
                                              ;
                       SEOR
                                              ; RETURN
   STATE28: MOV
                       STATE, #37
                                             ;STATE: =37
            YOM
                       ABORT,#36
                                             ;ABORT:=36
            CLR
                       MV05
                                             CLR VOS MASK
            CLR
                       MYC5
                                             CLR VC5 MASK
            SETB
                       CENO
                                             ; PRESS. CONTROL (ON)
            MOV
                       STPO, PATH
                                             SET ATM. SETPOINT
            CLR
                       FO
                                             CLEAR HOLD FLAG
            LJMP
                       SEOR
                                             :RETURN
  STATE29: MOV
                      C,LSC1
                                            :TEST V1 CLOSED-
           ANL
                       C,LSC2
                                            ; AND V2 CLOSED-
           ANL
                      C, LSC3
                                            ;AND V3 CLOSED-
           ANL
                      C,LSC4
                                            :AND V4 CLOSED-
           ANL
                      C, LSC5
                                            ; AND V5 CLOSED-
           ANL
                      C, LSC6
                                            ; AND V6 CLOSED-
           ANL
                      C, LS07
                                            ; AND V7 OPEN-
           ANL
                      C,LSC8
                                            ;AND V8 CLOSED-
           ANL
                      C,SUC2
                                            ; AND SU2 PUSHED
           JNC
                      S291
                                            : IF RESET
           YOK
                      STATE, #2
                                               STATE:=2
           MOV
                      ABORT, #0
                                               ABORT:=0
           YOM
                      STAT, $00H
                                              RESET STATUS
          MOV
                     CC00, #00H
                                              RESET ALARM LIGHTS
          MOV
                     CC01,#01H
                                              RESET RUN LIGHTS
          CLR
                     FO
                                              CLEAR HOLD FLAG
          SJMP
                     S292
S291:
                                              END
          YOM
                     CTRL, #00H
                                           ; ELSE, RESET CONTROLS
          MOV
                     TCEN. #00H
                                              RESET ALARMS
          HOV
                     MSKO, #OOH
                                              RESET CLOSED MASKS
          YOU
                     MSX1, #00H
                                              RESET OPEN MASKS
          YOK
                     MSK2, #00H
                                              RESET HISC. MASKS
          YOU
                                           ;
                     CC02, #40H
          HOV
                                           ÷
                                              RESET VALVES
                     CC03, $00H
                                              RESET MISC. OUTPUTS
                                          ;
         MOV
                     CC01, #00H
                                              TURN CYCLE LIGHTS OFF
         SETB
                                          ;
                    LT02
                                              EVAC-FAIL(ON)
         SETB
                    FO
S292:
                                              SET HOLD FLAG
         NOP
         LJMP
                                              END
                    SEOR
                                          :RETURN
```

	•		· ·
STATES	o: nov	C, LSC1	:
	ANL		:TEST V1 CLOSED
	ANL	C,LSC2	; AND V2 CLOSED-
		C, LSC3	; AND V3 CLOSED-
	ANL.	C, LSC4	;AND V4 CLOSED-
•	ANL	C, LSC5	; AND V5 CLOSED-
•	ANL	C,LSC6	· AND VS CLOSED- (
	ANL	C,LS07	; AND V6 CLOSED-
	ANL	C, LSC8	; AND V7 OPEN-
	ANL	C.SUC2	:AND V8 CLOSED-
•	JNC	S301	AND SU2 PUSHED
	MOV		; IF RESET
	MOV	STATE,#2	; STATE:=2
•	nov	ABORT, #0	; ABORT:=0
	MOV	STAT, #00H	; RESET STATUS
		CCOO, #00H	RESET ALARM LIGHTS
	MOV	CC01,#01H	RESET RUN LIGHTS
	CLR	FO	; CLEAR HOLD FLAG
	SJMP	S302	; END
S301:	MOV	CTRL, #OOH	
	MOV	TCEN, #OOH	ELSE, RESET CONTROLS
	MOV	MSKO, #OOH	; RESET ALARMS
•	MOV	MSK1, 400H	; RESET CLOSED MASKS
	MOV	MSK2, #00H	. ; RESET OPEN MASKS
•	MOV	CC03 140H	; RESET MISC. MASKS
	MOV	CC02, #40H	; RESET VALVES
	MOV	CC03, \$00H	; RESET MISC. OUTPUTS
	SETB	CC01, #00H	: TURN CYCLE LIGHTS OFF
		LT03	; FILL-FAIL(ON)
S302:	SETB	FO	; SET HOLD FLAG
2304:	NOP	•	; END
	LJMP	SEQR	RETURN
		•-) A S I U K M

```
STATE31: HOV
                         C,LSC1
                                               :TEST V1 CLOSED-
              ANL
                         C.LSO2
                                               ;AND V2 OPEN-
              anl
                         C,LSC3
                                               ; AND V3 CLOSED
              ANL
                         C.LSC4
                                               :AND V4 CLOSED-
              ANL
                         C. LSC6
                                               ;AND V6 CLOSED-
              ANL
                         C,LSC7
                                              .; AND V7 CLOSED-
              ANL
                         C,LSC8
                                             . ; AND V8 CLOSED-
              ANL
                         C.SUC2
                                               AND SUZ PUSHED
              JNC
                         S311
                                             FIF RESET
              HOV
                         STATE, #20
                                             .. ;. STATE:=20
             YOU
                        ABORT,#31
                                                 ABORT:=32
                                              :
             MOV
                        STAT, #00H
                                                 RESET STATUS
                                              :
             MOV
                        MSKO, #OEFH
                                                 SET ALL CLOSED MASKS
                                              •
             HOV
                        MSK1, #OEFH
                                              ;
                                                 SET ALL OPEN MASKS
             MOV
                        MSK2, #001H
                                                 SET MISC. MASKS
                                              ;
             MOV
                        CC00, 400H
                                                 RESET ALARM LIGHTS
             MOV
                        CC01,#22H
                                             . . RESET RUN LIGHTS
             CLR
                        FO
                                                 CLEAR HOLD FLAG
             SJMP
                        S312
   S311:
                                                 END
             YOM
                        CTRL, $03H
                                              ; ELSE, RESET CONTROLS
             HOV
                        TCEN, #00H
                                                 RESET ALARMS
             MOV
                        MSKO, #00H
                                                RESET CLOSED MASKS
RESET OPEN MASKS
             MOV
                        MSK1,#00H
                                              ;
             MOV
                        MSK2, #00H
                                                RESET MISC. MASKS
                                             ï
             MOV
                        CC02, #02H
                                                RESET ALL VALVES
                                             į
             YOM
                       CC03, #01H
                                                RESET MISC. OUTPUTS
                                             ;
             SETR
                       LTDA
                                                STERIL-FAIL(ON)
                                             ;
            SETB
                       FO
  S312:
                                             :
                                                SET HOLD FLAG
            NOP
                                                END
            LJMP
                       SEOR
                                             RETURN
  STATE32: NOV
                       C.LSC5
                                             TEST V5 CLOSED
            ANL
                       C,SUC2
                                             :AND SU2
            JNC
                       S321
                                             ; IF (V5 CLOSED & SU2 PUSHED)
            MOV
                       STATE, #19
                                             ; STATE:=19
            MOV
                       ABORT,#32
                                                ABORT:=32
            CLR
                       FO
                                                CLEAR HOLD FLAG
            SJMP
                      S322
                                                END
  S321:
            YON
                      CC02, #02H
                                            ; ELSE, RESET ALL VALVES
           SETB
                      FO
                                               SET HOLD FLAG
  S322:
           NOP
                                            : END
           LJMP
                      SEQR
                                            ; RETURN
 STATE33: NOV
                      C,SUC2
                                            :TEST SU2
           JNC
                      S331
           MOV
                                            : IF PUSHED
                      STATE, #23
                                               STATE:=23
           HOV
                      ABORT, #33
                                               ABORT:=34
           CLR
                      FO
                                               CLEAR HOLD FLAG
           SJMP
                      S332
 S331:
                                            ; END .
           SETB
                      FO
 S332:
                                            ; ELSE, SET HOLD FLAG
           NOP
                                            ; END
           LJMP
                      SEQR
                                            ; RETURN
STATE34: MOV
                     C.SUC2
                                           :TEST SU2
          JNC
                     5341
                                           : IF PUSHED
          MOV
                     STATE, #25
                                              STATE:=25
          MOV
                     ABORT,#35
                                              ABORT:=35
          LCALL
                     DCTO
                                              DECREMENT PURGE COUNT
          CLR
                    FO
                                              CLEAR HOLD FLAG
          SJMP
                    5342
S341:
                                              END
         SETB
                   . F0
                                          ; ELSE, SET HOLD FLAG
S342:
         NOP
                                             END
         LJMP
                    SEOR
                                          :RETURN
```

```
STATE35: MOV
                    C.LSC1
                                         :TEST V1 CLOSED-
          ANL
                    C,LS02
                                         ; AND V2 OPEN-
          ANL
                    C.LSC3
                                         :AND V3 CLOSED-
          ANL
                    C,LSC4
                                         ;AND V4 CLOSED-
                    C,LSC5
          ANL
                                        AND V5 CLOSED-
          ANL
                    C,LSC6
                                        :AND V6 CLOSED-
                    C,LSC7
          ANL
                                         :AND V7 CLOSED-
                                        ;AND V8 CLOSED-
                    C,LSC8
          ANL
          JNC
                    S351
                                         ; IF RESET
          YOR
                                         ; STATE:=20
                    STATE, $20
          MOV .
                    ABORT,#32
                                         ; ABORT:=32
          MOV
                    STAT, $00H
                                        ; RESET STATUS .
          MOV
                    MSKO, #OEFH
                                         ; SET ALL CLOSED MASKS
                    MSK1, #OEFH
          YOM
                                         ; SET ALL OPEN MASKS
          MOV
                    MSK2, 4001H
                                         ;
                                            SET MISC. MASKS
          VOM
                    CC00, #00H
                                         ;
                                            RESET ALARM LIGHTS
                                            RESET RUN LIGHTS
          VOM
                    CCO1, #22H
                                         ;
          YOM
                    STPO, PSP1
                                            LOAD PRESS. SETPOINT
                                         .
                   CENO
          SETB
                                            ENABLE PRESSURE CONTROL
          VOM
                    MTMO, #PN2T
                                            SET PRESSURE TIMER
                                         ÷
          CLR
                    TFL4
                                            CLEAR TIMER FLAG
                                         ;
          CLR
                    FO
                                            CLEAR HOLD FLAG
                                         ;
          SJMP
                    S352
                                            END
 S351:
          MOV
                    CTRL, #03H
                                         ; ELSE, RESET CONTROLS
                    TCEN. 400H
          VOM
                                           RESET ALARMS
                                         ;
          VOM
                                            RESET CLOSED MASKS
                    MSKO, #00H
                                         ;
          MOV
                    MSK1, #00H
                                         ; RESET OPEN MASKS
          MOV
                                         ; RESET MISC. MASKS
                    MSK2, #00H
                                         ; RESET ALL VALVES
          MOV
                    CCO2, #02H
          MOV
                    CC03, #01H
                                            RESET MISC. OUTPUTS
                                         ;
          SETB
                    FO
                                            SET HOLD FLAG
 S352:
          NOP
                                            END
                                        . ;
          LJMP
                                         ; RETURN
                    SEQR
 STATE36: MOV
                    C.SUC2
                                        :TEST SU2
          JNC
                    S361
                                         ; IF PUSHED
          MOV
                                         ; STATE:=26
                    STATE, $26
          VOM
                    ABORT,#37
                                            ABORT:=37
          CLR
                    FO
                                            CLEAR HOLD FLAG
          SJMP
                    5362
                                         ; END
S361:
          SETB
                    FO
                                         ; ELSE, SET, HOLD FLAG
5362:
          NOP
                                         ; END
          LJMP
                    SEQR
                                         ; RETURN
```

```
STATE37: CLR
                                           CLEAR CARRY
          MOV
                     A, & PATM
                                           GET ATM SETPOINT
          SUBB
                     A, ADIO
                                           SUBTRACT PRESSURE
          JC
                     5371
                                           ; IF PRESS.GT.ATM
          MOV
                     STATE, #38
                                              STATE:=38
          HOY
                     ABORT, #0
                                              ABORT:=0
          CLR
                     MVC7
                                             CLEAR VC7 MASK
          CLR
                     MVO7
                                            CLEAR VOT MASK
          CLR
                     CENO
                                             N2 LOOP(OFF)
          CLR
                     CTRO
                                             N2 OUTP(OFF)
          CLR
                     VV05
                                             CLOSE N2 VALVE
          SETB
                     VV07
                                             OPEN V7
          CLR .
                     LT16
                                              PURGE-IN-PROGRESS(OFF)
          SETB
                     LT17
                                             REMOVE-LOAD(ON)
          CLR
                     FO
                                             CLEAR HOLD FLAG
          SJMP
                    5372
                                             END
S371:
          SETB
                     FQ
                                          ; ELSE, SET HOLD FLAG
S372:
         NOP
                                             END
         LJMP
                    SEOR
                                          :RETURN
STATE38: JNB
                    SUC2.5381
                                          ; IF SU2 PUSHED
         MOV
                    STATE, $0
                                             STATE:=0 (RESET)
         VOM
                    ABORT, 40
                                             ABORT:=0
         CLR
                    FO
                                             CLEAR HOLD FLAG
         SJMP
                    S382
                                             END
S381:
         SETB
                    F0
                                          ELSE, SET HOLD FLAG
5382:
         NOP
                                             END
         LJMP
                    SEQR
                                          ; RETURN
;
```

In the foregoing specification, the invention has been described with reference to a specific exemplary embodiment thereof. It will, however, be evident that various modifications and changes may be made thereunto without parting from the broader spirit and scope of the invention as set forth in the appended claims. For example, as will be appreciated by those of ordinary skill in the art familiar with this specification, the apparatus

disclosed herein may be suitable for use in connection with various types of gaseous treatment systems, such as those which employ toxic gases, e.g., without limitation, bleaching gases, fumigants, sterilants, etc. The specification and drawings are, accordingly, to be regarded in an illustrative rather than a restrictive sense.

-86-

WHAT IS CLAIMED IS:

1. Apparatus for treating articles with a gas comprising: chamber means for receiving an article to be treated;

means for supplying the gas to the chamber means comprising valve means coupled to the chamber means for supplying the gas to the chamber means, means for removing the gas from the chamber means after a predetermined time interval, electronic control means receiving a plurality of electrical signals associated with ones of measured parameters from said chamber means for controlling said valve means and said means for removing, said electronic control means comprising computer means for cycling said apparatus through a plurality of states in accordance with a predetermined sequence of instructions, said computer means including means for aborting the operation of said apparatus to one of a plurality of defined failure states in response to a failure of said apparatus, said selected failure state dependent on the state in said cycle in which the failure occurred.

- The apparatus recited in claim 1 wherein said gas is a sterilizing gas, whereby said article is sterilized by said gas.
- 3. The apparatus recited in claim 2, wherein said sterilizing gas is generated from at least two component parts, and further including first means for receiving a first component part of the gas, second means for receiving a second component part of the gas, means for allowing said first and second component parts to react with each other to generate said sterilizing gas, said means for allowing being controlled by said computer means in response to the measurement of selected ones of said plurality of measured parameters.

- 4. The apparatus recited in claim 3, further comprising valve means for supplying a relatively stable gas to said chamber means.
- 5. The apparatus recited in claim 3, further comprising valve means for supplying filtered air to said chamber means.
- 6. The apparatus recited in claim 3, further comprising valve means for supplying water vapor to said chamber means to affect the humidity level in said chamber.
- 7. The apparatus recited in claim 3 wherein said plurality of measured parameters include the temperature, pressure and humidity in said chamber means and the concentration of said sterilizing gas in said chamber means.
- 8. The apparatus recited in claim 2 wherein said sterilizing gas comprises chlorine dioxide.
- 9. The apparatus recited in claim 3 wherein said sterilizing gas comprises chlorine dioxide and said first component comprises chlorine gas and said second component comprises sodium chlorite.
- 10. The apparatus recited in claim 2 wherein said means for removing comprises vacuum pump means and additional valve means.
- 11. The apparatus recited in claim 2, wherein said valve means comprises first and second switch means, said first switch means indicating when said valve means is open and said second switch means indicating when said valve means is closed, said first and second switch means being in opposite states such that when said first switch means is closed, said second switch means is open.

- 12. The apparatus recited in claim 2, wherein said electronic control means comprises memory means, and further comprising means for receiving input signals from said valve means indicative of the closed or open condition of said valve means and means for transmitting output signals to said valve means to selectively open or close said valve means, image signals of said input and output signals being stored in said memory means.
- 13. The apparatus recited in claim 12, further comprising means for disabling said output signals from being transmitted to said valve means except when an enabling signal is issued by said computer means.
- 14. The apparatus recited in claim 2, further comprising means for monitoring for proper operation of said computer means, said monitoring means issuing a disabling signal to prevent actuation of said valve means in the event of a failure of said computer means.
- 15. The apparatus recited in claim 12, further comprising mask means stored in said memory means, said computer means comparing said image signals of said input and output signals and generating an alarm signal if said input and output image signals do not agree in response to the setting of a bit in said mask means.
- 16. The apparatus recited in claim 11, further comprising means for monitoring the state of said first and second switch means, and further comprising means for generating an alarm signal if said first and second switch means are not in the proper states.
- 17. The apparatus recited in claim 2, wherein said valve means moves between a first state and a second state in response to instructions from said computer means, and further comprising timer means for generating an alarm

signal if said valve means does not move from said first to second state in a predetermined time interval.

- 18. The apparatus recited in claim 2, further comprising means for cycling said apparatus to a further defined state once one of said defined failure states is reached.
- 19. Apparatus for treating articles with a gas comprising:
 first means for receiving a first component;
 second means for receiving a second component, said
 first and second components, when reacted together,
 forming said gas;

means for reacting said two components together for forming said gas;

first valve means for supplying said gas to said chamber means to treat said article in said chamber means;

means for removing said gas from said chamber means; electronic controller means for controlling said means for reacting, means for supplying and means for removing comprising computer means executing a predetermined sequence of steps so as to cycle said apparatus through a series of successive states defining a cycle in which said article is treated by said gas and wherein said gas is thereafter removed from said chamber means so as to render said chamber means within acceptable standards of safety.

- 20. The apparatus recited in claim 19 wherein said gas is a sterilizing gas, whereby said article is sterilized by said gas.
- 21. The apparatus recited in claim 20 wherein said computer means comprises means for receiving a plurality of electrical signals associated with ones of measured parameters from said chamber means for controlling the operation of said means for reacting, means for supplying

and means for removing.

- 22. The apparatus recited in claim 21, wherein said means for reacting comprises second valve means for allowing said first and second components to react with each other to generate said sterilizing gas, said second valve means being controlled by said computer means in response to the measurement of selected ones of said plurality of measured parameters.
- 23. The apparatus recited in claim 22, further comprising valve means for supplying a relatively stable gas to said chamber means.
- 24. The apparatus recited in claim 22, further comprising valve means for supplying filtered air to said chamber means.
- 25. The apparatus recited in claim 22, further comprising valve means for supplying water vapor to said chamber means to affect the humidity level in said chamber.
- 26. The apparatus recited in claim 22 wherein said plurality of measured parameters include the temperature, pressure and humidity in said chamber means and the concentration of said sterilizing gas in said chamber means.
- 27. The apparatus recited in claim 20 wherein said sterilizing gas comprises chlorine dioxide.
- 28. The apparatus recited in claim 27 wherein said sterilizing gas comprises chlorine dioxide and said first component comprises chlorine gas and said second component comprises sodium chlorite.
- 29. The apparatus recited in claim 20 wherein said means for removing comprises vacuum pump means and additional valve

means.

- 30. The apparatus recited in claim 20, wherein said valve means comprises first and second switch means, said first switch means indicating when said valve means is open and said second switch means indicating when said valve means is closed, said first and second switch means being in opposite states such that when said first switch means is closed, said second switch means is open.
- 31. The apparatus recited in claim 20, wherein said electronic control means comprises memory means, and further comprising means for receiving input signals from said valve means indicative of the closed or open condition of said valve means and means for transmitting output signals to said valve means to selectively open or close said valve means, images of said input and output signals being stored in said memory means.
- 32. The apparatus recited in claim 31, further comprising means for disabling said output signals from being transmitted to said valve means except when an enabling signal is issued by said computer means.
- 33. The apparatus recited in claim 20, further comprising means for monitoring for proper operation of said computer means, said monitoring means issuing a disabling signal to prevent actuation of said valve means in the event of a failure of said computer means.
- 34. The apparatus recited in claim 31, further comprising mask means stored in said memory means, said computer means comparing said images of said input and output signals and generating an alarm signal if said input and output images do not agree in response to the setting of a bit in said mask means.

- 35. The apparatus recited in claim 30, further comprising means for monitoring the state of said first and second switch means, and further comprising means for generating an alarm signal if said first and second switch means are not in the proper states.
- 36. The apparatus recited in claim 20, wherein said valve means moves between a first state and a second state in response to instructions from said computer means, and further comprising timer means for generating an alarm signal if said valve means does not move from said first to second state in a predetermined time interval.
- 37. The apparatus recited in claim 20, wherein said computer means includes means for aborting the operation of said apparatus to one of a plurality of defined failure states in response to a failure of said apparatus, said selected failure state dependent on the state in said cycle in which the failure occurred.
- 38. The apparatus recited in claim 37, further comprising means for resetting said apparatus to a further defined state once one of said defined failure states is reached.
- chamber means for receiving articles to be treated;
 means for supplying the gas to the chamber means
 comprising valve means coupled to the chamber means for
 supplying the gas to the chamber means, means for
 removing the gas from the chamber means after a
 predetermined time interval, electronic control means
 receiving a plurality of electrical signals associated
 with ones of measured parameters from said chamber means
 for controlling said valve means and said means for
 removing, said electronic control means comprising
 computer means for cycling said apparatus through a
 plurality of states in accordance with a predetermined

sequence of instructions,

said computer means including memory means, and further comprising means for receiving input signals from said valve means indicative of the closed or open condition of said valve means and means for transmitting output signals to said valve means to selectively open or close said valve means, image signals of said input and output signals being stored in said memory means,


mask means being stored in said memory means, said computer means comparing said image signals of said input and output signals and generating an alarm signal if said input and output image signals do not agree in response to the setting of a bit in said mask means.


- 40. The apparatus recited in claim 39, wherein said gas is a sterilizing gas, whereby said article is sterilized by said gas.
- 41. The apparatus recited in claim 40, wherein said sterilizing gas is generated from at least two component parts, and further including first means for receiving a first component part of the gas, second means for receiving a second component part of the gas, means for allowing said first and second component parts to react with each other to generate said sterilizing gas, said means for allowing being controlled by said computer means in response to the measurement of selected ones of said plurality of measured parameters.
- 42. The apparatus recited in claim 41, further comprising valve means for supplying a relatively stable gas to said chamber means.
- 43. The apparatus recited in claim 41, further comprising valve means for supplying filtered air to said chamber means.


- 44. The apparatus recited in claim 41, further comprising valve means for supplying water vapor to said chamber means to affect the humidity level in said chamber.
- 45. The apparatus recited in claim 41, wherein said plurality of measured parameters include the temperature, pressure and humidity in said chamber means and the concentration of said sterilizing gas in said chamber means.
- 46. The apparatus recited in claim 40 wherein said sterilizing gas comprises chlorine dioxide.
- 47. The apparatus recited in claim 41 wherein said sterilizing gas comprises chlorine dioxide and said first component comprises chlorine gas and said second component comprises sodium chlorite.
- 48. The apparatus recited in claim 40 wherein said means for removing comprises vacuum pump means and additional valve means.
- 49. The apparatus recited in claim 40, wherein said valve means comprises first and second switch means, said first switch means indicating when said valve means is open and said second switch means indicating when said valve means is closed, said first and second switch means being in opposite states such that when said first switch means is closed, said second switch means is open.
- 50. The apparatus recited in claim 40, further comprising means for disabling said output signals from being transmitted to said valve means except when an enabling signal is issued by said computer means.
- 51. The apparatus recited in claim 40, further comprising means for monitoring for proper operation of said computer means, said monitoring means issuing a disabling

signal to prevent actuation of said valve means in the event of a failure of said computer means.

- 52. The apparatus recited in claim 49, further comprising means for monitoring the state of said first and second switch means, and further comprising means for generating an alarm signal if said first and second switch means are not in the proper states.
- 53. The apparatus recited in claim 40, wherein said valve means moves between a first state and a second state in response to instructions from said computer means, and further comprising timer means for generating an alarm signal if said valve means does not move from said first to second state in a predetermined time interval.
- 54. The apparatus recited in claim 40, wherein said computer means further comprises means for aborting the operation of said apparatus to one of a plurality of defined failure states in response to a failure of said apparatus, said selected failure state dependent on the state in said cycle in which the failure occurred.
- 55. The apparatus recited in claim 39, further comprising means for cycling said apparatus to a further defined state once one of said defined failure states is reached.







4/25

· 1	1				
ADDRESS	DESCRIPTION	A	A	a AL	3 ^A 12
ØØ-FF	INTERNAL RAM	_			
0000-0FFF	INTERNAL ROM	Ø	Ø	ø	ø
1000-1FFF	EXTERNAL ROM	Ø	Ø	ø	/
2000-203F	EXTERNAL SRAM	Ø	Ø	/	Ø
4000-400F	CLOCK	Ø	/	ø	ø
6000-6007	A/D READ	Ø	/	/	Ø
COOO	X00-X07)	1	7	Ø	Ø
CØ01	X10-X17	/	/	Ø	Ø
CØ02.	X20-X27 DIN	1	/	Ø	Ø
C003	x30-x37)	/	/	ø	ø
Εφφφ	Y00-Y07)	1	/	/	φ
E001	Y10-Y17 (/	1	/	Ø
E002	Y20 Y27 DOU	/	1	1	ø
EØØ3	Y30 Y37	1	1	1	Ø
E Ø Ø 4	WATCHDOG - RESET	,	/	/	Ø

FIG. 3A (a)


FIG.	FIG.
3A(a)	3A(b)

FIG.3A

ADDRESS BUS BITS

**************************************	1/3	
A ₁₁ A ₁₀ A ₉ A ₈	A7 A6 A5 A4	A3 A2 A1 A0
	A7 A6 A5 A4	A3 A2 A1 A0
AII AIOA9 A8	A7 A6 A5 A4	A3 A2 A1 A0
AII AIO A9 A8	A7 A6 A5 A4	A3 A2 A1 A6
Ø Ø Ø Ø	Ø Ø A5 A4	A3 A2 A1 A6
<i>• • • • • •</i>	Ø Ø Ø Ø	A3 A2 A1 A6
0000	0 0 0 0	Ø A2 A1 A6
ø ·	0	0000
Ø Ø	Ø	0001
Ø	Ø Ø	Ø Ø 1 Ø Ø Ø 1 1
Ø	Ø	0000
ø	Ø	0001
ø	ø	0010
	Ø	0011
P	Ø	0100

FIG. 3A (b)




7/25


OLTI - DOOR OPEN
OLT2 - EVAC FAIL
OLT3 - FILL FAIL
OLT4 - STERIL FAIL
OLT5 - PURGE FAIL
OLT6 - LOAD UNSTERILE
OLT7
OLT8

O LTII-READY FOR CYCLE O LTI2-CYCLE IN PROGRESS O LTI3-EVAC IN PROGRESS O LTI4-FILL IN PROGRESS O LTI5-STERIL IN PROGRESS O LTI6-PURGE IN PROGRESS O LTI7-REMOVE LOAD O LTI8-

START CYCLE

ABORT-RESET





10/25

FIG. 7


F/G. 7


				l				29					1	30)	ı			P	RO 31	CE:	SS	•
				0	1	2	3	4	5	6	7	8	9			12	13	14	15	_	17	18	19
01591445	DOOR OPEN EVAC FAIL FILL FAIL STERIL FAIL PURGE FAIL LOAD MASTER	LTØ1 LTØ3 LTØ4 LTØ5 LTØ6 LTØ7 LTØ8		1000000	10000000	0000000	0000000	0000000	0000000	0000000	0000000	0000000	0000000	0000000	0000000	0000000	0000000	0000000	0000000	0000000	00000000	0000000	0000000
5	READY CYCLE EVAC FILL STERIL PURGE ROM LOAD	LT14	001	000000	000000	000000	100000	100000	100000	110000	110000	110000	101000	101000	101000	100100	100100	100700	100100	100100	01000100	100010	100010
V A L V E S	MAIN VAC GAS ENABLE VAC CTRL GAS CTRL N2 CTRL H2O CTRL ATM VENT VAC VLV.	VVØ1 VVØ3 VVØ4 VVØ5 VVØ6 VVØ7 VVØ8	CC02	000001	000001	000001	000000	000000	000000	000000	000000	000000	00000000	0000000	1000000	10000	100000	100000	10000	10000	01100000	11000	10000
		HTOI	6	0	0	0	0	/	C	C (۱ء	C	C	C	C	C (C	C	00		


F/G. 7α


STATES 32 33 3435 20 21 22 23 24 25 2	627 28 20 20 21 22 23 34 35
000000000000000000000000000000000000000	000 00000000000000000000000000000000000
0 0 0 0 0 0 0	00 0000000000
000000000000000000000000000000000000000	00 000000000000000000000000000000000000
00000000	00 0000000000
000000000000000000000000000000000000000	10 0011101000 10 0000100000 10 000000000
00000100	0 0011111000


FIG. 7b























DATA MEMORY MAP

FIG. 12A





FIG. 12B









FIG. 15







	CONTACT OUTPUT STATUS]	(WCO)	1	CONTACT OUTPUTS
28:	CCOØ	1 .	WRITE	E000:	YØ
29:	CCOI		CONTACT OUTPUTS	EDD1:	×1
2A:	CCO2	1	0017013	E002:	Y2
28:	CCO3	1		E003:	Y3

FIG. 20






FIG. 21



SUBCTITUTE SHEET






FIG. 24

SUBSTITUTE SHEET

INTERNATIONAL SEARCH REPORT

			International Application No PCT/U	S86/00258		
		ON OF SUBJECT MATTER (if several classi				
		ational Patent Classification (IPC) or to both Nati	lonal Classification and IPC			
-		305D 7/06 12/111		,		
II. FIELD						
		Minimum Documer	ntation Searched 4			
Classificati	ion System		Classification Symbols			
U.S. 422/3, 27, 28, 29, 34, 37, 110, 111, 114, 116, 295, 305; 364/413, 499, 500						
		Documentation Searched other to the Extent that such Documents	han Minimum Documentation are included in the Fields Searched •			
		CONSIDERED TO BE RELEVANT 14				
ategory *	l C	tation of Document, 16 with Indication, where app	ropriate, of the relevant passages 17	Relevant to Claim No. 15		
$\frac{X}{Y}$		A, 4,067,691, (McGady et al) See the entire document.	4,067,691, (McGady et al), 10 January 1978, ee the entire document.			
Y,P		A, 4,504,442, (Rosenblatt et See column 3, lines 27-47; c and column 5, lines 5-46.	3-10, 19-38, 41-48			
Y		7, 21- ee column 2, lines 35-38 and 57-69 and column 1 lines 1-13.				
Y		A, 4,431,159, (Stubbs), 14 1 column 4, lines 4-7.	11-13, 15-17 30-32, 34-36 39-55			
Y		12, 13, 2-36. (Grudowski), 13 September 1983, 12, 13, 31, 32, 39-55				
¥ .		A, 3,910,761, (Hopkins), 07 column 10, lines 1-5.	October 1975, See	18, 38, 55		
"A" doccor "E" ear fill "L" docwh cit: "O" docott	cument di nsidered i riler document watch is cit- cument watch is cit- ation or o cument re- cument per cument per than the	iblished prior to the international filing date but e priority date claimed	"I" later document published after to or priority date and not in conflicted to understand the principle invention "X" document of particular relevant cannot be considered novel or involve an inventive step "Y" document of particular relevant cannot be considered to involve document is combined with one ments, such combination being in the art. "&" document member of the same p	ct with the application but a critheory underlying the ce; the claimed invention cannot be considered to ce; the claimed invention an inventive stop when the or more other such documberous to a person skilled		
Date of th	·	ON Completion of the International Search ³	Date of Mailing of this International Se			
22	April	1986	0 6 MAY 1986			
		ning Authority 1	Signature of Authorized Officer 20	/		
ISA/	US		Brion P. Heaney			

III. DOCU	CUMENTS CONSIDERED TO BE RELEVANT (CONTINUED FROM THE SECOND SHEET)					
Category •	Citation of Document, 16 with indication, where appropriate, of the relevant passages 17	Relevant to Claim No 18				
A	GB, A, 2,052,800, (Cook et al.), 28 January 1981, See entire document.	1–55				
A	US, A, 4,164,538, (Young et al.), 14 August 1979, See entire document.	1–55				
A	US, A, 4,239,731, (Gillis et al.), 16 December 1980, See entire document.	1-55				
A	US, A, 4,261,950, (Bainbridge et al.), 14 April 1981, See entire document.	1-55				
A	US, A, 4,294,804, (Baran), 13 October 1981, See entire document.	1-55				
A	US; A, 4,372,916, (Chamberlain et al.), 08 February 1983, See entire document.	1-55				
A	US, A, 4,447,399, (Runnells et al.), 08 May 1984, See entire document.	1-55				
A	US, A, 4,457,892, (Young), 03 July 1984, See entire document.	1-55				
1	dicroprocessors & Microsystems, Volum 3, No. 8, published October 1979 (Great Britain), R.N. Mewis, "Triplicated Microprocessor Controlled Automatic Shutdown System", see pages 347 to 351.	1-55				
	; :	. •				
ı						